EN

ES

FR

DE

IT

PT

SE

NL

DK

PL

CZ

HU

# PJX Series Balances

## User Guide

## 1. SAFETY INFORMATION

This manual contains installation, operation and maintenance instructions for PJX Series balance. Please read this manual completely before installation and operation.

### Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

|                       |                                                                                                                                            |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>WARNING</b>        | For a hazardous situation with medium risk, possibly resulting in injuries or death if not avoided.                                        |
| <b>CAUTION</b>        | For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or injuries if not avoided. |
| <b>Attention Note</b> | For important information about the product<br>For useful information about the product                                                    |

### Warning Symbols

General Hazard

Electrical Shock Hazard

Alternating current

Direct current

### Safety Precautions

**CAUTION:** Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

- Verify that the AC adapter's input voltage range and plug type are compatible with the local AC mains power supply.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- Do not position the balance such that it is difficult to reach the power connection.
- The balance is for indoor use only. Do not operate the equipment in hazardous or unstable environments.
- Operate the equipment only under ambient conditions specified in these instructions.
- Do not drop loads on the pan.
- Use the balance only in dry locations.
- Disconnect the equipment from the power supply when cleaning.
- Use only approved accessories and peripherals.
- Service should only be performed by authorized personnel.

### Intended Use

Use the instrument exclusively for weighing as described in the operating instructions. Any other type of use and operation beyond the limits of technical specifications without written consent from OHAUS, is considered as not intended.





This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use.

If the instrument is not used according to these operating instructions, the intended protection of the instrument may be compromised and OHAUS assumes no liability.

## 2. INSTALLATION

### 2.1 Select the Location

Avoid heat sources, rapid temperature changes, air current or excessive vibrations. Allow sufficient space.


### 2.2 Leveling

Be sure the balance is level before it is used or after its location is changed.

The PJX balance has a level bubble in a small round window beside the display.

To level the balance, adjust the 4 leveling feet until the bubble is centered in the circle.

Please refer to the right figure for leveling.

### 2.3 Connecting Power and Acclimatising the Balance

Connect the DC output connector to the power receptacle on the rear of the balance. Then connect the AC adapter plug to a suitable electrical outlet.

#### Acclimatising

It is suggested that the balance should not be used until it has been connected to power and acclimatised to the environment for a certain period of time. In the case of a balance with the precision above 0.1 mg, the acclimatisation time should be 1.5 hours; in the case of balance with the precision of 0.01 mg, the acclimatisation time should be more than 4 hours.

### 2.4 Connecting the Interface

The PJX balance has two connectivities, RS232 port and Mini DIN.

Use the RS-232 port to connect either to a computer or a printer with a standard (straight-through) serial cable.

Mini DIN is used to connect with another accessory, as an alternative option for RS232.

Interface connections on the rear of the balance.

Mini DIN

RS232

RS232: Used to connect to PC or Printer

**Note:** See the Printing section for Connecting, Configuring and Testing the Printer / Computer Interface.


## 2.5 Initial Calibration

For best results, the instrument's calibration should be checked prior to first use. If adjustment is needed, refer to the Calibration section of the instruction manual.

## 3. OPERATION

### 3.1 Overview of Display, Home Screen

#### CONTROLS

#### CONTROL FUNCTIONS

| Button                                                                                                                                                   | <div>Zero </div> <div>Yes</div>                                          | <div>Print Unit</div> <div>No</div>                                                                                                                           | <div>Function Mode</div> <div>Back</div>                                                                                                     | <div>Tare Menu-Cal</div> <div>Exit</div>                                                                                                                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Primary Function (Short Press)</b>                                                                                                                    | <b>On / Zero</b> <ul style="list-style-type: none"> <li>If the balance is Off, turns on the balance.</li> <li>If the balance is On, sets zero.</li> </ul> | <b>Print</b> <ul style="list-style-type: none"> <li>Sends the current displayed value to the serial interface.</li> </ul>                                     | <b>Function</b> <ul style="list-style-type: none"> <li>Operation is dependent on the application mode.</li> </ul>                            | <b>Tare</b> <ul style="list-style-type: none"> <li>Performs tare operation.</li> </ul>                                                                           |
| <b>Secondary Function (Long Press)</b><br><br><b>(Extended Press)</b> | <b>Off</b> <ul style="list-style-type: none"> <li>Zeros the current value.</li> <li>None</li> </ul>                                                       | <b>Unit</b> <ul style="list-style-type: none"> <li>Changes weighing units.</li> <li>None</li> </ul>                                                           | <b>Mode</b> <ul style="list-style-type: none"> <li>Changes application mode.</li> <li>None</li> </ul>                                        | <b>Menu-Cal</b> <ul style="list-style-type: none"> <li>Enters the main menu. Calibration is the first sub-menu.</li> <li>Views the preset Tare value.</li> </ul> |
| <b>Menu Function (Short Press)</b>                                                                                                                       | <b>Yes</b> <ul style="list-style-type: none"> <li>Accepts the current (blinking) setting on the display.</li> </ul>                                       | <b>No</b> <ul style="list-style-type: none"> <li>Rejects the current (blinking) setting on the display.</li> <li>Increments a value being entered.</li> </ul> | <b>Back</b> <ul style="list-style-type: none"> <li>Reverts back to previous menu item.</li> <li>Decrements a value being entered.</li> </ul> | <b>Exit</b> <ul style="list-style-type: none"> <li>Immediately exits the sub-menu.</li> <li>Aborts a calibration in progress.</li> </ul>                         |

## DISPLAY

### 3.2 Principal Functions

**Zeroing:** Press **Zero** to zero the balance.

**Taring:** With no load on the pan, press **Zero** to set the display to zero. Place an empty container on the pan and press **Tare**. Add objects to the container and its net weight is displayed. After the container and the objects are removed, the tare value will be displayed as a negative number. Press **Tare** to clear.

**weighing:** First press **Zero** to set the display to zero. Place an object on the pan. The display indicates the gross weight.

## 4. MAINTENANCE

### 4.1 Cleaning

**WARNING:** Disconnect the balance from the power supply before cleaning. Make sure that no liquid enters the interior of the balance.

Clean the Balance at regular intervals.

Housing surfaces may be cleaned with a lint-free cloth slightly dampened with water or a mild cleaning agent.

Glass surfaces may be cleaned with a commercial glass cleaner.

**Attention:** Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

### 4.2 Troubleshooting

If the troubleshooting section does not resolve your problem, contact an Authorized Ohaus Service Agent. Please visit our website [www.ohaus.com](http://www.ohaus.com) to locate the Ohaus office nearest you.

## 5. TECHNICAL DATA

The technical data is valid under the following conditions:

- Indoor use only
- Altitude: Up to 2000 m
- Specified Temperature range: 10°C to 30°C
- Humidity: maximum relative humidity 80% for temperatures up to 30°C, decreasing linearly to 50% relative humidity at 40°C
- Mains supply voltage fluctuations: up to  $\pm 10\%$  of the nominal voltage
- Installation category II
- Pollution degree: 2

| Model                         | PJX203M                                                                    | PJX803M | PJX1203M | PJX1602M | PJX2502M |
|-------------------------------|----------------------------------------------------------------------------|---------|----------|----------|----------|
| Capacity (ct)                 | 200                                                                        | 800     | 1200     | 1600     | 2500     |
| Capacity (g)                  | 40                                                                         | 160     | 240      | 320      | 500      |
| Readability d (ct)            | 0.001                                                                      | 0.001   | 0.001    | 0.01     | 0.01     |
| Readability d (g)             | 0.0001                                                                     | 0.0001  | 0.0001   | 0.001    | 0.001    |
| Verification Interval e (g)   | 0.001                                                                      | 0.001   | 0.001    | 0.01     | 0.01     |
| Class                         | II                                                                         | I       | I        | II       | II       |
| Repeatability (std. dev.) (g) | 0.0001                                                                     | 0.0001  | 0.0001   | 0.001    | 0.001    |
| Linearity (g)                 | 0.0002                                                                     | 0.0002  | 0.0002   | 0.002    | 0.002    |
| Power Supply                  | Power input: 100-240V ~ 200mA 50-60Hz 12-18VA<br>Power output: 12 VDC 0.5A |         |          |          |          |

| Model                         | PJX822M                                                                    | PJX1602M | PJX2202M | PJX4202M | PJX5202M | PJX4201M |
|-------------------------------|----------------------------------------------------------------------------|----------|----------|----------|----------|----------|
| Capacity (g)                  | 820                                                                        | 1600     | 2200     | 4200     | 5200     | 4200     |
| Readability d (g)             | 0.01                                                                       | 0.01     | 0.01     | 0.01     | 0.01     | 0.1      |
| Verification Interval e (g)   | 0.1                                                                        | 0.1      | 0.1      | 0.1      | 0.1      | 0.1      |
| Class                         | II                                                                         | II       | II       | II       | II       | II       |
| Repeatability (std. dev.) (g) | 0.01                                                                       | 0.01     | 0.01     | 0.01     | 0.01     | 0.1      |
| Linearity (g)                 | 0.02                                                                       | 0.02     | 0.02     | 0.02     | 0.02     | 0.2      |
| Power Supply                  | Power input: 100-240V ~ 200mA 50-60Hz 12-18VA<br>Power output: 12 VDC 0.5A |          |          |          |          |          |

| Model                         | PJX203<br>PJX203/E                                                         | PJX803<br>PJX803/E | PJX1203<br>PJX1203/E | PJX1603<br>PJX1603/E | PJX2503<br>PJX2503/E |
|-------------------------------|----------------------------------------------------------------------------|--------------------|----------------------|----------------------|----------------------|
| Capacity (ct)                 | 200                                                                        | 800                | 1200                 | 1600                 | 2500                 |
| Capacity (g)                  | 40                                                                         | 160                | 240                  | 320                  | 500                  |
| Readability d (ct)            | 0.001                                                                      | 0.001              | 0.001                | 0.005                | 0.005                |
| Readability d (g)             | 0.0001                                                                     | 0.0001             | 0.0001               | 0.001                | 0.001                |
| Repeatability (std. dev.) (g) | 0.0001                                                                     | 0.0001             | 0.0001               | 0.001                | 0.001                |
| Linearity (g)                 | 0.0002                                                                     | 0.0002             | 0.0002               | 0.002                | 0.002                |
| Power Supply                  | Power input: 100-240V ~ 200mA 50-60Hz 12-18VA<br>Power output: 12 VDC 0.5A |                    |                      |                      |                      |

| Model                         | PJX822<br>PJX822/E                                                         | PJX1602<br>PJX1602/E | PJX2202<br>PJX2202/E | PJX3202 |
|-------------------------------|----------------------------------------------------------------------------|----------------------|----------------------|---------|
| Capacity (g)                  | 820                                                                        | 1600                 | 2200                 | 3200    |
| Readability d (g)             | 0.01                                                                       | 0.01                 | 0.01                 | 0.01    |
| Repeatability (std. dev.) (g) | 0.01                                                                       | 0.01                 | 0.01                 | 0.01    |
| Linearity (g)                 | 0.02                                                                       | 0.02                 | 0.02                 | 0.02    |
| Power Supply                  | Power input: 100-240V ~ 200mA 50-60Hz 12-18VA<br>Power output: 12 VDC 0.5A |                      |                      |         |

| Model                         | PJX4202<br>PJX4202/E                                                       | PJX5202<br>PJX5202/E | PJX2201<br>PJX2201/E | PJX4201<br>PJX4201/E | PJX6201<br>PJX6201/E |
|-------------------------------|----------------------------------------------------------------------------|----------------------|----------------------|----------------------|----------------------|
| Capacity (g)                  | 4200                                                                       | 5200                 | 2200                 | 4200                 | 6200                 |
| Readability d (g)             | 0.01                                                                       | 0.01                 | 0.1                  | 0.1                  | 0.1                  |
| Repeatability (std. dev.) (g) | 0.01                                                                       | 0.01                 | 0.1                  | 0.1                  | 0.1                  |
| Linearity (g)                 | 0.02                                                                       | 0.02                 | 0.2                  | 0.2                  | 0.2                  |
| Power Supply                  | Power input: 100-240V ~ 200mA 50-60Hz 12-18VA<br>Power output: 12 VDC 0.5A |                      |                      |                      |                      |

6. COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

| Mark                                                                             | Standard                                                                                                                                                                                                                |
|----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | This product complies with the EU Directives 2014/30/EU (EMC), 2014/35/EU (LVD) and 2014/31/EU (NAWI). The EU Declaration of Conformity is available online at <a href="http://www.ohaus.com/ce">www.ohaus.com/ce</a> . |
|  | EN 61326-1, AS/NZS 61010-1                                                                                                                                                                                              |
|  | CAN/CSA-C22.2 No. 61010-1<br>UL Std. No. 61010-1                                                                                                                                                                        |

Verified weighing instruments

When the instrument is used in trade or a legally controlled application it must be set up, verified and sealed in accordance with local weights and measures regulations. It is the responsibility of the purchaser to ensure that all pertinent legal requirements are met.

Weighing Instruments verified at the place of manufacture bear the following supplementary metrology marking on the descriptive plate.

Weighing Instruments to be verified in two stages have no supplementary metrology marking on the descriptive plate. The second stage of conformity assessment must be carried out by the applicable weights and measures authorities.

If national regulations limit the validity period of the verification, the user of the weighing instrument must strictly observe the re-verification period and inform the weights and measures authorities.

As verification requirements vary by jurisdiction, the purchaser should contact their local weights and measures office if they are not familiar with the requirements.

|                                                                                     |                                                                                                                                                                                                               |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. |
|  | For disposal instructions in Europe, refer to <a href="http://www.ohaus.com/weee">www.ohaus.com/weee</a>                                                                                                      |

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P/N 30431477 C © 2018 Ohaus Corporation, all rights reserved.
Ohaus Corporation