


Ziegelei 1 D-72336 Balingen

Phone: +49-[0]7433-9933-0 Fax: +49-[0]7433-9933-149 E-Mail: info@kern-sohn.com Internet: www.kern-sohn.com

Operating manual Multi-function weighing balance





MWA-M / MWA-PM-BA-e-1614



KERN MWA 300K1M / MWA 300K1PM Version 1.4 04/2016 Operating manual Multi-function weighing balance

Contents

1	Technical data	4
2 2.1	Declaration of conformity Explanation of the graphic symbols	
3	Appliance overview	7
4	Keyboard overview	9
5	Overview of display	10
6 6.1 6.2 6.3 6.4 6.5	Basic instructions Specific function Proper use Improper Use Warranty Monitoring of Test Resources	11 11 12 12
7 7.1 7.2 7.3	Basic Safety Precautions Pay attention to the instructions in the Operation Manual Personnel training Preventing contamination	13 13
8 8.1 8.2 8.3 8.3.1 8.4	Electromagnetic compatibility (EMC) General hints Electromagnetic interferences Electromagnetic noise immunity Crucial features of performance Minimum distances	14 15 16 18
9 9.1 9.2	Transport and storage Testing upon acceptance Packaging / return transport	19
10 10.1 10.2 10.3 10.4 10.4.1 10.5 10.6 10.7	Unpacking, Setup and Commissioning Installation Site, Location of Use Unpacking Scope of delivery Balance assembly and installation Application examples: Mains connection Battery operation is possible by obtaining an optional battery power pack Battery operation	20 21 21 21 21 22 23 24
10.8	Initial Commissioning	26

11	Operation	
11.1	Weighing	
11.1.1	Weighing with wheelchair	
11.2	Taring	27
11.3	HOLD function	
11.4	Display additional last decimal place	29
11.5	Calculation of the Body Mass Index	29
11.5.1	Calculating Body Mass Index	
11.5.2	Classification of BMI values	
11.6	Automatic switch-off function "AUTO OFF"	
11.7	Display background illumination	
12	Menu	
12.1	Navigation in the menu	
12.2	Menu overview	
13		
13.1	Data output RS 232	
13.1	Pin allocation of balance output bushing:	
13.2 13.3	Technical data	
13.5	Printer operation	
14	Error messages	37
15	Servicing, maintenance, disposal	
15.1	Cleaning	
15.2	Cleaning / disinfecting	
15.3	Servicing, maintenance	
15.4	Disposal	
16	Instant help	39
17	Verification	40
17.1	Verification validity period (current status in G)	
18	Adjustment	
-	•	

1 Technical data

KERN	MWA 300K1M	MWA 300K1PM	
Display	6-c	ligit	
Weighing range (max)	300) kg	
Minimum load (Min)	2	kg	
Verification value (e)	10	0 g	
Reproducibility	0.1	kg	
Linearity ±	0.1	kg	
Display	LCD with 25n	nm high digits	
Recommended adjustment weight, (Class)		300 kg (M1)	
Stabilization time (typical)	2 s	ec.	
Warm-up time	10 min		
Operating temperature	- 10° C + 40° C		
Humidity of air	max. 80 % (not condensing)		
Electric Supply	Input voltage 100V-240V AC, 50 / 60 Hz		
Auto Off	After 3, 5, 15, 30 min witho	ut load change (adjustable)	
Dimensions fully assembled (D x W X H) mm	1190 x 1140 x 80	1190 x 1140 x 1150	
Weighing plate (W x D) mm	880 x 840		
Weight kg (net)	72	76	
Rechargeable battery operation	Optional		
Batteries	6 x 1.5	5 V AA	
Data interface provided as standard	RS 232 C		


2 Declaration of conformity

To view the current EC/EU Declaration of Conformity go to:

<u>www.kern-sohn.com/ce</u>

• The scope of delivery for verified weighing balances (= conformityrated weighing balances) includes a Declaration of Conformity. Solely these weighing balances are classified as medical devices.

2.1 Explanation of the graphic symbols





"Observe operating instructions"

"Observe operating instructions"



Identification of manufacturer of medical product including address

Kern & Sohn GmbH D–72336 Baligen,Germany www.kern-sohn.com



"Electro-medical appliance" with attachment for type B




Device protection category II



Dispose of old appliances separately from your household waste!

Instead, take them to communal collection points.



Temperature limit indicating the upper and the lower limit (storage temperature on packaging) (Temperature serving as an example)



Display of supply voltage for scales with polarity display.



Safety instruction: Do not place hands beneath the ramp.

3 Appliance overview






Weighing balance suitable as wheelchair or stand assist balance.


- 1. Display Unit
- Handrail optional MWA-A01with tripod (cannot be retrofitted) MWA-A02 without tripod
- 3. Ramps
- 4. Bubble level
- 5. Carrying handle





English

4 Keyboard overview



Button	Designation	Function
ON OFF	ON/OFF button	Turn on/off
HOLD PRINT G	Hold-key Print-key	Hold function /determination of stable weighing value Print (press and hold for a while): Data transmission via interface RS 232
BMI K	BMI key	Determination of body mass index In the menu: • Confirm selection For numeric entry: • Confirm numerical value
F	Function key	In menu: Call up menu Select menu items For numeric entry: Increase numerical value
→0← ←	Zeroing key	Weighing scale will be reset to "0.0" For numeric entry: • Change decimal place
TARE	TARE key	Tare balance

5 Overview of display

Display	Designation	Description
	Stability display	Scales are in a steady state
→0 ←	Zeroing display	Should the balance not display exactly zero despite empty weighing plate, press the button. Your balance will be set to zero after a short standby time.
NET	Net weight display	Illuminated when net weight is displayed Illuminated after weighing scale was tared
GROSS	Gross weight display	Illuminated when gross weight is displayed
HOLD	HOLD function	HOLD function active
BMI	BMI function	Illuminated while BMI function is enabled
	Rechargeable /Battery display	Shows the capacity of the rechargeable / batteries

6 Basic instructions



Weighing instruments have to be verified for the purposes stated below in accordance with Directive 2009/23/EC. Article 1, paragraph 4. "Determination of mass in the practice of medicine that is, weighing patients for reasons of medical supervision during medical surveillance, examination and treatment."

6.1 Specific function

- **Indication** Determining the body weight in the medical practice area.
 - Use as "non-standalone weighing scale", that is, a person steps carefully onto the weighing platform's centre. Once a steady display value is shown, you can read the weight value.
 - For weighing with a wheelchair, the wheelchair and the person sitting on it are pushed over the ramp onto the centre of the weighing platform or, in the case of an electric wheelchair, the patient and wheelchair move independently onto the weighing platform. Once a steady display value is shown, you can read the weight value.
- **Contraindication** No contraindication known

6.2 Proper use

This weighing balance is designed for determining the weight of persons in a standing or sitting position in medical treatment rooms. The balance is suitable for recognising, preventing and controlling illnesses.

On multi-function balances, the weighed person should step onto the centre of the weighing platform and remain standing without moving. For weighing with a wheelchair, the wheelchair should be placed in the centre of the weighing platform.

As soon as a stable weighing value is reached the weighing value can be read. The weighing scale is designed for continuous duty.



The weighing platform may only be stepped on by persons capable of standing on both feet on the weighing platform.

The balance should be checked for correct condition prior to each utilisation by a person familiar with proper operation of the balance.

6.3 Improper Use

Do not use these scales for dynamic weighing processes.

Do not leave permanent load on the weighing pan. This may damage the measuring system.

Impacts and overloading exceeding the stated maximum load (max) of the weighing plate, minus a possibly existing tare load, must be strictly avoided. This could cause damage to the balance.

Never operate balance in explosive environment. The serial version is not explosion protected. It should be noted that a flammable mixture of anaesthetics and oxygen or laughing gas may occur.

The structure of the balance may not be modified. This may lead to incorrect weighing results, safety-related faults and destruction of the balance.

The balance may only be used according to the described conditions. Other areas of use must be released by KERN in writing.

6.4 Warranty

Warranty claims shall be voided in case

- Our conditions in the operation manual are ignored
- The appliance is used outside the described uses
- The appliance is modified or opened
- Mechanical damage and damage caused by media, liquids,
- Natural wear and tear
- The appliance is improperly set up or incorrectly electrically connected
- The measuring system is overloaded
- Dropping the balance

6.5 Monitoring of Test Resources

In the framework of quality assurance the measuring-related weighing properties of the balance and, if applicable, the testing weight, must be checked regularly. The responsible user must define a suitable interval as well as type and scope of this test. Information is available on KERN's home page (<u>www.kern-sohn.com</u> with regard to the monitoring of balance test substances and the test weights required for this. In KERN's accredited DKD calibration laboratory test weights and balances may be calibrated (return to the national standard) fast and at moderate cost.

12

7 Basic Safety Precautions

7.1 Pay attention to the instructions in the Operation Manual



Carefully read this operation manual before setup and commissioning, even if you are already familiar with KERN balances.

7.2 Personnel training

The medical staff must apply and follow the operating instructions for proper use and care of the product.

7.3 Preventing contamination

The prevention of cross-contamination (fungal skin infections,.....) requires regular cleaning of the weighing platform. Recommendation: after a weighing procedure that could potentially result in contamination (e. g. after weighing that involves direct skin contact).

8 Electromagnetic compatibility (EMC)

8.1 General hints



The installation and use of the electrical multi-function balances MWA requires special precautionary measures as outlined in the EMC information below.

This device complies with the limits set for medical electrical devices of group 1, class B (as per EN 60601-1-2).

Electromagnetic compatibility (EMC) describes a device's ability to perform reliably within an electromagnetic environment without causing inadmissible electromagnetic interference at the same time. Amongst other things, such disturbances may be emitted by connecting cables or the air.

Inadmissible disturbances from the environment may result in incorrect displays, inaccurate measured values or incorrect behaviour of the multi-function balances MWA. The multi-function balances MWA, too, may in certain cases cause such malfunctions in other devices. To eliminate problems of that kind, we recommend you to take one or several of the measures listed below:

- Change the alignment or distance of the device to the source of failure.
- Install and use the multi-function balances MWA in one and the same place.
- Connect the multi-function balances MWA to a separate power source.
- For further questions please contact our customer services.

Disturbances may be caused by improper modification or add-ons to the device or not recommended accessories (such as power units or connecting cables). The manufacturer will not be responsible for these. Modifications may also result in a loss of authorisation relating to the use of the device.



Devices emitting high frequency signals (mobile telephones, radio transmitters, radio receivers) may cause interference in the multi-function balances MWA. Do not use these near multi-function balances MWA. Chapter 8.4 contains details about recommended minimum distances.

8.2 Electromagnetic interferences

Guidelines and manufacturer's declaration – electromagnetic interferences

The multi-function balances MWA are intended for operation in the kind of electromagnetic environment stated below. The customer or user of the multi-function balances MWA should ensure that they are operated in this kind of environment.

Emitted interference measurements	Conformity	Electromagnetic environment - guideline
HF emissions as per CISPR 11 / EN 55011	Assembly 1	The multi-function balances MWA apply HF energy exclusively for their internal function. Their HF emission therefore is very low and it is highly unlike to interfere with adjacent electronic devices.
HF emissions as per CISPR 11 / EN 55011	Class B	The multi-function balances MWA are designed for use in all equipment including those in living
Emission of harmonics acc. to IEC 61000-3-2	Class A	areas and those connected directly to the public power grid that also supplies buildings used for living purposes.
Emission of voltage fluctuations / flicker	Conforms with	
acc. to IEC 61000-3-3		

Do not stack or operate multi-function balances MWA immediately next to or together with other devices. If this type of operation is necessary, observe the multi-function balances MWA to ensure normal operation in such an arrangement.

8.3 Electromagnetic noise immunity

Guidelines and manufacturer's declaration - electromagnetic noise immunity

The multi-function balances MWA are intended for operation in the kind of electromagnetic environment stated below. The customer or user of the multi-function balances MWA should ensure that they are operated in this kind of environment.

Noise immunity tests	IEC 60601 test level	Conformity	Electromagnetic environment - guideline
Discharge static electricity (DSE) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV ± 8 kV	Floors should be made of wood or concrete or tiled with ceramic tiles. If floors are covered with synthetic material, relative air humidity must be at least 30%.
Fast transient electrical disturbances / bursts acc. to IEC 61000-4-4	± 2 kV for power lines <u>+</u> 1 kV for input and output lines	± 2 kV <u>+</u> 1 kV	The quality of the supply voltage should match that of the typical business or hospital environment.
Impulse voltages / surges acc. to IEC 61000-4-5	± 1 kV voltage Live wire - live wire ± 2 kV voltage Live wire - earth	± 1 kV Inapplicable	The quality of the supply voltage should match that of the typical business or hospital environment.
Voltage dips, short-term disruptions and fluctuations in supply voltage acc. to IEC 61000-4-11	< 5 % U _T (> 95 % dip of U _T) for ½ period 40 % U _T (> 60 % dip of U _T) for 5 periods 70 % U _T (> 30 % dip of U _T) for 25 periods < 5 % U _T (> 95 % dip of U _T) for 5 s	Compliance with requirements under all postulated conditions Controlled switch off Return to undisturbed situation after user intervention.	The quality of the supply voltage should match that of the typical business or hospital environment. If the user demands continued functioning of multi-function balances MWA even after the occurrence of a power failure, we recommend using uninterruptible means of power supply or batteries for the multi-function balances MWA.
Magnetic field for supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m 50/60 Hz	Magnetic fields for the supply frequency should match the typical values found in the particular business or hospital environment.
NOTE U_T equals A	C line voltage prior to a	oplication of te	est level.

Guidelines and manufacturer's declaration - electromagnetic noise immunity

The multi-function balances MWA are intended for operation in the kind of electromagnetic environment stated below. The customer or user of the multi-function balances MWA should ensure that they are operated in this kind of environment.

Noise immunity tests	IEC 60601 test level	Conformity	Electromagnetic environment - guideline
Conducted HF disturbance variables acc. to IEC 61000-4- 6	3 V _{rms} 150 kHz to 80 MHz	3 V	Do not use portable or mobile radio sets nearer to the multi-function balances MWA or their wires than the distance recommended as safety distance which is calculated according to the equation relevant for its transmission frequency.
Emitted HF disturbance variables	3 V _{rms} 80 MHz to 2.5 GHz	3 V/m	Recommended safety distance: $d = 1.2\sqrt{P}$
Acc. to IEC 61000-4-			$d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz
3			$d = 2.3\sqrt{P}$ for 800 MHz up to 2.5 GHz
			Use P as rated capacity of radio transmitter in Watt (W) as per details given by the radio transmitter manufacturer and d as recommended safety distance in metres (m).
		$((\bullet))$	The field intensity of stationary radio transmitters should for all frequencies be lower according to an in situ ^a examination than the conformity level. ^b
			Interference may occur near devices bearing the symbol below.
NOTE 2 T	ligher frequency range a hese guidelines may not he spread of electromag eflections in buildings, ob	t be applicable ir netic variables i	n all cases. s influenced by absorption and
mobile radio sets, a reliably predicted in stationary transmitt If the measured fiel conformity level ab normal operation. It	amateur radio stations, A n advance. To determine ers, you should consider Id intensity at the location ove, you should observe f you observe unusual fe	M and FM radio the electromag r a study of elect n where the bala the multifunctio eatures of perform	ase stations of wireless telephones and o and television stations cannot be netic environment in respect of tromagnetic phenomena at the location. ance is to be used exceeds the onal balances MWA in order to ensure mance you may have to take additional tion for the multifunctional balances
^b For a frequency rar	nge of 150 kHz to 80 MH	Iz field intensity	should be less than 3 V/m.

8.3.1 Crucial features of performance

Note:



The multifunctional balances MWA do not have any crucial features of performance as per IEC 60601-1. The system may be subject to interference by other devices even if these devices conform to current emission requirements as per CISPR.

8.4 Minimum distances

Recommended safety distance between portable and mobile HF telecommunication devices and multi-function balances MWA

The multi-function balances MWA are designed for use in an electromagnetic environment in which HF disturbance variables are controlled. The customer or user of the multi-function balances MWA can help avoiding electromagnetic disturbances by keeping the minimum distance between portable and mobile HF telecommunication devices (transmitters) and the multi-function balances MWA – depending on the output performance of the communication device, as stated below.

Rated capacity of transmitter W	The safety distance depends on the transmission frequency m		
	150 kHz up to 80 MHz $d = 1.2\sqrt{P}$	80 MHz up to 800 MHz $d = 1.2\sqrt{P}$	800 MHz up to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.00

For transmitters with a maximum rated capacity not stated in the table above you can calculate the recommended safety distance in metres (m) yourself by using the equation belonging to each column, whereby P equals the maximum rated capacity of the transmitter in Watt (W) as per details provided by the transmitter manufacturer.

NOTE 1Higher frequency range applies to 80 MHz and 800 MHz.NOTE 2These guidelines may not be applicable in all cases.The spread of electromagnetic variables is influenced by absorption and
reflections in buildings, objects and humans.

9 Transport and storage

9.1 Testing upon acceptance

When receiving the appliance, please check packaging immediately, and the appliance itself when unpacking for possible visible damage.

9.2 Packaging / return transport



- ⇒ Keep all parts of the original packaging for a possibly required return.
- ⇒ Only use original packaging for returning.
- ⇒ Prior to dispatch disconnect all cables and remove loose/mobile parts.
- ⇒ Reattach possibly supplied transport securing devices.
- Secure all parts such as the weighing platform, power unit etc. against shifting and damage.

10 Unpacking, Setup and Commissioning

10.1 Installation Site, Location of Use

The balances are designed in a way that reliable weighing results are achieved in common conditions of use.

You will work accurately and fast, if you select the right location for your balance.

On the installation site observe the following:

- Place scales on a stable, even surface;
- Avoid extreme heat as well as temperature fluctuation caused by installing next to a radiator or in the direct sunlight;
- Protect the balance against direct draughts due to open windows and doors;
- Avoid jarring during weighing;
- Protect the balance against high humidity, vapours and dust;
- Do not expose the device to extreme dampness for longer periods of time. Non-permitted condensation (condensation of air humidity on the appliance) may occur if a cold appliance is taken to a considerably warmer environment. In this case, acclimatize the disconnected appliance for ca. 2 hours at room temperature.
- Avoid static charge of the balance and of the person to be weighed.
- Avoid contact with water.

Major display deviations (incorrect weighing results) may be experienced should electromagnetic fields (e.g. due to mobile phones or radio equipment), static electricity accumulations or instable power supply occur. Change location or remove source of interference.


10.2 Unpacking

Remove the individual components of the balance or the complete balance from the packaging with care and install at the intended location. When using the power pack, ensure that the power cable does not produce a risk of stumbling.


10.3 Scope of delivery

Serial accessories:

- Balance
- Mains adapter (in conformity with EN 60601-1)
- Operating manual
- Wall bracket
- Small components



10.4 Balance assembly and installation




- ⇒ Level balance with foot screws until the air bubble of the water balance is in the prescribed circle.
- ⇒ Check levelling regularly.

The weighing balance to be used as wheelchair balance is delivered fully assembled. 1 handrail with tripod (MWA-A01) or 1 handrail without tripod (MWAS-A02) are optionally available for additional use as stand assist balance (See chap.3 device overview). The handrails are also suitable for use in connection with wheelchairs.

10.4.1 Application examples:

Weighing balance used as			
Wheelchair balance	stand assist balance with handrail MWS- A01		

Assembly handrail MWA-A02 (optional):





Note:

After fitting the attachments and before using the weighing balance, check all attachments for firm fit.

10.5 Mains connection

Power is supplied by the external power unit which also serves to isolate the mains supply from the scale. The stated voltage value must be the same as the local voltage.

Always use genuine approved KERN power pack units as per EN 60601-1 directive.

The small sticker attached to the side of the display unit indicates the power port:

_	(F	
	1	-	-	
	ł	Ļ	,	

The LED remains illuminated as long as the weighing scale remains connected to the mains.

The LED display informs you during loading about the loading status of the rechargeable battery.

- green: Rechargeable battery completely reloaded
- blue: Charging storage battery

English



10.6 Battery operation is possible by obtaining an optional battery power pack

Open the battery compartment cover (1) at the base of the display unit and insert the rechargeable battery. Charge the battery for at least 12 hours before initial use.

The appearance of the symbol in the weight display indicates that the battery is almost exhausted. The weighing scale will remain ready for operation for a few more minutes before switching off in order to save battery. Load rechargeable battery.



Voltage has dropped below prescribed minimum.



Rechargeable battery very low.



Rechargeable battery completely reloaded

Right underneath the display there is a LED with the symbol **1**. If the LED lights green, the rechargeable battery is fully charged. If it is lighting blue, it is being charged.

If the balance is not used for a longer time, take out the rechargeable battery and store it separately. Leaking liquid could damage the balance.

10.7 Battery operation

As an alternative to rechargeable battery operation, the balance may also be operated with 6x AA batteries.

Open battery compartment cover (1) at the lower side of the display unit and insert batteries according to the example below. Lock the battery cover again. If the

batteries are empty, in the balance display appears the symbol batteries. To save battery power, the balance switches off automatically (see chap.11.6 Auto off).



Capacity of batteries exhausted.



Batteries will soon be flat.



Batteries are completely charged

Insert batteries

Remove battery compartment cover	
Connect battery holder to housing contact acc. to illustration	
Insert battery holder	0 + 4 + 1 +
Insert batteries into battery compartment and lock with battery compartment cover.	

10.8 Initial Commissioning

In order to obtain exact results with the electronic balances, your balance must have reached the operating temperature (see warming up time chap. During this warming up time the balance must be connected to the power supply (mains, accumulator or battery) and be switched on.

The accuracy of the balance depends on the local acceleration of gravity. The value of gravity acceleration is shown on the type plate.

11 Operation

11.1 Weighing

ZERO	00
GROSS	

⇒

Start balance by pressing
The balance will carry out a self-test
The scales are ready for operation as soon as the weight
display for "0.0 kg" has appeared.

ON

	 However, you can reset the weighing scale to zero by pressing
İ	the key.
	,

- ⇒ Have person stand in the centre of the scales. Wait until the standstill display "STABLE" appears, then read the weighing result.
- If the person is heavier than the weighing range, "OL" (=overload) will appear in the display.

- \Rightarrow Place wheelchair and patient in the centre of the weighing balance.
- ⇒ Engage brakes on wheelchair




Do not leave the patient unattended

- \Rightarrow Wait until patient is sitting still and read weighing value 1.
- ➡ Loosen the brakes and carefully pull off the transportation stretcher/wheelchair with patient
- ⇒ Then weigh the wheelchair without the patient and deduct the result from weighing value 1, from this follows the weight for the patient.


11.2 Taring

The tare weight of any preloads can be deducted by pressing a button so that the actual weight of the person is displayed in subsequent weighings.



11.3 HOLD function

The balance has an integrated standstill function (mean value calculation). With this function it is possible to weigh people accurately even if they do not stand still on the weighing plate.



11.4 Display additional last decimal place

(short-time, additional decimal place)


Press and hold for about 2 s whilst weighed result is being shown. The second decimal place will be shown for approx. 5 s.

This value is not considered as verified and may not be used for the purpose of a verified balance.

11.5 Calculation of the Body Mass Index

You need to know a person's body height before you can calculate the BMI for that person. This should be known.

11.5.1 Calculating Body Mass Index



English

- Reliable calculation of BMI is restricted to a body height of 100 cm to 200 cm and a weight of >10 kg.
- If weighing has to take place under unsteady conditions, you can be stabilise the display by applying the Hold function.


11.5.2 Classification of BMI values

Weight classification for adults over 18 years of age using the BMI in accordance with WHO, 2000 EK IV and WHO 2004.

Categorie	BMI (kg/m²)	Risk of diseases associated with overweight
Underweight	< 18.5	low
Normal weight	18.5 – 24.9	Average
Overweight	<u>></u> 25.0	
Pre-adipose	25.0 – 29.9	A bit high
Adipose degree I	30.0 - 34.9	Increased
Adipose degree II	35.0 - 39.9	High
Adipose degree III	<u>≥</u> 40	Very high

11.6 Automatic switch-off function "AUTO OFF"


The weighing scale will switch off automatically after the allotted time as long as neither the display unit nor the weighing plate is operated.



0.0 kg

GROSS

11.7 Display background illumination



12 Menu



Access to service menu "tCH" is locked in verified balances. To disable the access lock, destroy the seal and actuate the adjustment switch. For position of adjustment switch, see chap. 17. **Attention**:

After destruction of the seal the weighing system must be re-verified by an authorised agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.

12.1 Navigation in the menu

Call up menu	⇒ In weighing mode, press and the first function [F1 oFF] will be displayed.
Select function	⇒ With help of , the individual functions can be selected one after the other.
Change settings	 ⇒ Confirm selected function by Logical. ⇒ Select desired setting by Logical and confirm with Logical setting, the balance returns to the menu.
Exit menu/ Return to weighing mode	⇒ Press , the balance will return to weighing mode.

12.2 Menu overview

Function	Settings	Description
Automatic cutout	oFF 0*	Automatic shutdown off
	oFF 3	Automatic shutdown after 3 min.
	oFF 5	Automatic shutdown after 5 min.
	oFF 15	Automatic shutdown after 15 min.
	oFF 30	Automatic shutdown after 30 min.
F 2 5 J E	oFF*	Not documented
Background	Prt	
of display	Pr ACC	
	bl on	Back lighting for display on
F3 6F	bl oFF	Display background illumination off
Background illumination of display	bl AU*	Backlighting for display will come on automatically as soon as the weighing scale is operated.
FYP-E Interface parameter		e by , then confirm with . Weight will be added to summation memory and printed
	P Prt	after pressing PRINT
	P Cont	Continuous data output
	Series	Not documented
	ASK	Remote control instructions: W: Send all weighing details S: Send stable weight value T: Taring Z: Zeroing
	P cnt 2	Not documented
	P Stab	Automatic data output of stable weighing values
	P Auto	Weighed result will be added automatically to summation memory and issued
	2. Baud rate The currently set ba	aud rate (b xxx) will be shown after the RS-232 mode was
	вмі	esired Baudrate by pressing and confirm by pressing and rate: 600, 1200, 2400, 4800, 9600

	3 Data	output for	mat		
	3. Data output format (Setting P Prt, P Auto, P Cont only!)				
	The currently set data output format will be displayed after the Baud rate has				
	been confirmed. Select desired format by and confirm with				
		Prt 0-3	Data output format, see chap	o. 13	
	only at setting P Prt, P				
		Cont 1	Default	Sd0 – on/off Continuous data output, selectable "Sending 0" yes / no	
	Only when set P Cont	Cont 2	Not documented		
		Cont 3	Not documented		
	 4. Printer type The currently set printer type will be displayed after has been confirmed. Select desired printer type by pressing and 		вмі		
	LP 50	·	Not documented		
	tPUP		Use this setting		
E [H Service menu	Pin		Password Input: press		
	Operate	e adjustme	nt switch; for position see cha	p.15	
P I SPJ	15*		-		
Display speed	30 60		Not documented		
	7.5		-		
	1.0				
	Adjustme	nt, see cha	p. 18		
PS [8L					
	·				
	tri*		Not documented		
P3 Pro	CoUnt		Not documented		
	rESEt		Reset weighing scale to fact	ory setting	
	SEtGrA		Not documented		

13 Data output RS 232

You can print weighing data automatically via the RS 232 interface or manually by HOLD PRINT O


via the interface according to the setting in the menu. pressing

This data exchange is asynchronous using ASCII - Code.

The following conditions must be met to provide successful communication between the weighing balance and the printer.

- Use a suitable cable to connect the weighing balance to the interface of the printer. Faultless operation requires an adequate KERN interface cable.
- Communication parameters (baud rate, bits and parity) of balance and printer must match. For detailed description of interface parameters see chap. 13.2)

13.1 Pin allocation of balance output bushing:



Pin 2: TXD - output
Pin 3: RXD - input
Pin 5: GND - signal ground

13.2 Technical data

Connection	9 pin d-subminiature bushing		
	Pin 2 output		
	Pin 3 input		
	Pin 5 signal earth		
Baud rate	Optional 600/1200/2400/4800/9600		
Parity	8 bits		

13.3 Printer operation

Printout examples:		
Prt		
0/2	60.0kg	
1/3	60.0kg 170.0cm 20.7BMI	

14 Error messages

Display

isplay	Description
Erry	Zero range exceeded
	 (on start-up or when pressing the key) Load on weighing pan
	 Excess load, during zero setting of weighing scale Incorrect adjusting process Fault on load cell
8008	Value outside the A/D converter rangeDamaged weighing cell

Damaged electronics •

Should other error messages occur, switch balance off and then on again. If the error message remains inform manufacturer.

15 Servicing, maintenance, disposal

15.1 Cleaning



Before any maintenance, cleaning and repair work disconnect the appliance from the operating voltage.

15.2 Cleaning / disinfecting

Clean weighing platform (such as seat) as well as casing with household detergents or commercially available disinfectants. Please follow manufacturer's instructions.

Do not use abrasive or aggressive cleaners such as spirits or alcohol or similar as they might damage the high-quality surface.

The prevention of cross-contamination (fungal skin infections,.....) requires regular cleaning of the weighing platform. Recommendation: after a weighing procedure that could potentially result in contamination (e. g. after weighing that involves direct skin contact).



Do not spray disinfectants onto appliance.

Make sure that disinfectant does not penetrate the interior of the appliance.

Remove dirt immediately.

15.3 Servicing, maintenance

The appliance may only be opened by trained service technicians who are authorized by KERN.

Disconnect the scales before opening.

15.4 Disposal

Disposal of packaging and appliance must be carried out by operator according to valid national or regional law of the location where the appliance is used.

16 Instant help

In case of an error in the program process, briefly turn off the balance and disconnect from power supply. The weighing process must then be restarted from the beginning.

Fault	Possible cause
The displayed weight does not glow.	• The balance is not switched on.
	 The mains supply connection has been interrupted (mains cable not plugged in/faulty).
	Power supply interrupted.
	 Rechargeable battery inserted incorrectly or empty
	No rechargeable battery inserted
The displayed weight is permanently changing	 Draught/air movement Table/floor vibrations
	 Table/fillor vibrations The weighing plate is in contact with foreign bodies or is not correctly positioned.
	 Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)
5 5	The display of the balance is not at zero
obviously incorrect	Adjustment is no longer correct.
	Great fluctuations in temperature.
	Warm-up time was ignored.
	 Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)
Should other error messages	occur, switch balance off and then on again. If the error

Should other error messages occur, switch balance off and then on again. If the error message remains inform manufacturer.

17 Verification

General introduction:

According to EU directive 2009/23/EC balances must be officially verified if they are used as follows (legally controlled area):

- a) For commercial transactions if the price of goods is determined by weighing.
- b) For the production of medicines in pharmacies as well as for analyses in the medical and pharmaceutical laboratory.
- c) For official purposes
- d) For manufacturing final packages

In cases of doubt, please contact your local trade in standard.

Verification notes:

An EU type approval exists for balances described in their technical data as verifiable. If a balance is used where obligation to verify exists as described above, it must be verified and re-verified at regular intervals.

Re-verification of a balance is carried out according to the respective national regulations. For verification validity period, s. chap. 17.1.

The legal regulation of the country where the balance is used must be observed!




Verification of the balance is invalid without the seal.

The seal marks attached on balances with type approval point out that the balance may only be opened and serviced by trained and authorised specialist staff. If the seal mark is destroyed, verification looses its validity. Please observe all national laws and legal regulations. In Germany a reverification will be necessary.

Balances with obligation to verify must be taken out of operation if:

- The weighing result of the balance is outside the error limit. Therefore, in regular intervals load balance with known test weight (ca. 1/3 of the max. load) and compare with displayed value.
- The reverification deadline has been exceeded.

Position adjustment switch and seals



- 1. Self-destroying seal mark
- 2. Cover
- 3. Adjustment switch

17.1 Verification validity period (current status in G)

Multi-function balances (including chair and wheelchair scales) in hospitals	4 years
Multi-function balances, when not located in hospitals (for example, doctor's offices and nursing homes)	unlimited
Baby weighing scales and mechanical birth weight scales	4 years
Bed scales	2 years
Scales in dialysis stations	unlimited

Rehab clinics and health authorities are treated as hospitals (4 years of verification validity)

Not treated as hospitals (verification validity not limited) are dialysis stations, nursing homes and doctor's surgeries.

(Data source : "Bureau of Standards News, Weighing Instruments in Medicine"}

English

18 Adjustment

As the acceleration value due to gravity is not the same at every location on earth, each display unit with connected weighing plate must be coordinated - in compliance with the underlying physical weighing principle - to the existing acceleration due to gravity at its place of location (only if the weighing system has not already been adjusted to the location in the factory). This adjustment process must be carried out for the first commissioning, after each change of location as well as in case of fluctuating environment temperature. To receive accurate measuring values it is also recommended to adjust the display unit periodically in weighing operation.

1	•	Prepare the required adjustment weight. The adjustment weight to be applied depends on the capacity of a weighing scale, see chap. 1. Carry out adjustment as closely as possible to admissible maximum load of weighing scale. Info about test weights can be found on the Internet at: http://www.kern-sohn.com.
	•	Observe stable environmental conditions. For warm-up time required for stabilisation see chap. 1.

\land	Access to service menu "tCH" is locked in verified balances. To disable the access lock, destroy the seal and actuate the adjustment switch. Position of the adjustment switch see chap. 17.
	Attention: After destruction of the seal the weighing system must be re-verified by an authorised agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.

Procedure:

