**Scout™**  
**USB Device Interface**  
INSTRUCTION MANUAL

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## INTRODUCTION

This Interface Kit is for use with OHAUS Scout STX, SPX, SKX and SJX series products.

The OHAUS USB Interface Kit is a unique solution to the problem of connecting a balance to a computer using a Universal Serial Bus (USB). The USB Interface kit includes a CD with the software drivers to create the required virtual port on the computer.

## SYSTEM REQUIREMENTS



- PC running Windows 98®, Windows 98SE®, Windows ME®, Windows 2000®, Windows XP® or Windows 7®
- Be forward compatible with iMac 10.7 (contact with OHAUS)
- Available USB port

## INTERFACE INSTALLATION

Install the Interface module on the mini DIN port (round connector) located at the rear side of the balance as shown. During this process, please make sure the balance is powered off.

Notes:

- Please make sure the small triangle mark (▲) on the round connector is at the bottom and then plug in.
- The appearance of your model may be different.

## SETUP

- A) The OHAUS Interface is preset to communicate using the following settings: **9600 baud, 8 bit, no parity, no handshake**. Determine the USB Device parameters required for the computer that is to be connected. If the parameters do not match, it will be necessary to change either the balance settings, or the computer settings.
- B) Upon installation of the interface, the balance will recognize the USB Device Interface and add relevant items to the menu. Configure the balance to the desired USB and printing parameters; refer to the balance Instruction Manual for assistance in using the menus.

### For SPX, SKX and SJX Balances:

Upon installation, the balance will show [ᄀᄀᄀᄀᄀ] during the power sequence.

<b>PRINT</b>	<b>USB.Dev</b>
Reset: <b>no</b> , yes	Reset: <b>no</b> , yes
Stable: <b>off</b> , on	Baud Rate: 1200... <b>9600</b> ...115200
A.Print: <b>off</b>	Parity: 7 even, 7 odd, 7 none, <b>8 none</b>
On.Stable interval	Stop bit: <b>1, 2</b>
accept	Handshake: <b>none</b> , XOn-XOff
continuous	Alternate command: Print (A...P...Z) Tare (A...T...Z) Zero (A...Z)
Content: Result (-> <b>off</b> , on)	End USB.Dev: Exit menu
Gross (-> <b>off</b> , on)	
Net (-> <b>off</b> , on)	
Tare (-> <b>off</b> , on)	
Header (-> <b>off</b> , on)	
Footer (-> <b>off</b> , on)	
Mode (-> <b>off</b> , on)	
Unit (-> <b>off</b> , on)	
Info (-> <b>off</b> , on)	
Layout: Format (->Single, <b>Multi</b> )	
Feed (-> <b>Line</b> , 4 Lines, Form)	
End Print: Exit menu	




**Note:** Items with bold font are default settings.

### SETUP DEFINITIONS

<b>Menu Settings</b>	<b>Explanation</b>
PRINT / A. Print – Off	Only sends data when PRINT is pressed.
PRINT / A.Print / On.Stable – Load, Load.Zero	Balance only transmits stable data. <i>Load:</i> Print stable value excluding zero <i>Load.Zero:</i> Print stable value including zero
PRINT / A.Print / interval - (xx) sec	Balance sends data every (xx) seconds.
PRINT / Auto Print – accept	Balance only sends stable accept data in Check mode.
PRINT / A.Print – Continuous	Balance repeatedly sends data as fast as possible.
PRINT / Content	On: print relevant content after weighing data Gross – on: G; Net – on: N; Tare – on: T; Info – on: balance will print application mode setting parameters
PRINT / Layout / Format - (Single, Multi)	<i>Single:</i> print all data in one line <i>Multi:</i> print all data in multiple lines
PRINT / Layout / Feed - (Line, 4 Lines, Form)	<i>Line:</i> feed one line after printing <i>4 Lines:</i> feed four lines after printing <i>Form:</i> feed one page after printing (move to the top of next page after printing)

**For STX Balances:**

Upon installation, the STX balance will show the USB icon () in the upper right corner.

Menu added	Function
Communication  --USB Device  --Baud Rate	Set the baud rate (bits per second). 
--Transmission	Set the data bits, stop bit, and parity. 
--Handshake	Set the flow control method. 
--Print Settings	For more info, please refer to the STX instruction manual. 

## VIRTUAL PORT SOFTWARE INSTALLATION

1. Insert the supplied CD into the CD drive.

Different versions of Windows® have slightly different steps to load the driver that is on the CD. In all versions the New Hardware Wizard guides you through the required steps to select the driver that is located on the CD.

Example of Windows XP Hardware Wizard

You may also get the latest drivers from this link:  
<http://www.ftdichip.com/Drivers/D2XX.htm>

2. After clicking Finish, the virtual port should be ready for use.

Windows® typically adds the virtual port in sequence after the highest number COM port. For example, on PC's equipped with up to 4 COM ports, the virtual port will be COM5.

When using the USB interface with programs that limit the number of COM port designations, it may be necessary to assign one of these port numbers to the new virtual port. This can be done in the Port Settings of the Device Manager utility, found in the Windows Control Panel.

## USB DEVICE OUTPUT

With the interface installed the balance will operate in several ways according to the settings.

Print format can be switched by xFMT user command.

**New Scout Print Format (Default Format):**

Output String (Non Check Weighing Applications):

[weight]	11 characters (right justified)
[space]	1 character
[unit]	5 characters (right justified)
[space]	1 character
[stability indicator]	1 character; "?" when unstable, space when stable
[space]	1 character
[T/N/G/PT]	2 characters (right justified)
[Term]	2 characters

**Note:** All of the fields have fixed length.

Output example:

```

*****192.21_****g_*_**
*****0.01_****g?_**

*****95.0_****g_*_N
*****169.6_****g_*_G
*****95.0_****g_*_N
*****74.6_****g_*_T
    
```

A.Print: off; Stable: off  
 Stable reading  
 unstable reading

A.Print: off; Stable: on  
 Content / Result -> on  
 Content / Gross -> on  
 Content / Net -> on  
 Content / Tare -> on

Output String (Check Weighing Application):

[weight]	11 characters (right justified)
[space]	1 character
[unit]	5 characters (right justified)
[space]	1 character
[stability indicator]	1 character; "?" when unstable, space when stable
[space]	1 character
[T/N/G/PT]	2 characters (right justified)
[space]	1 character
[application status]	6 characters (right justified)
[Term]	2 characters

Example:

```

*****192.21_****g_*_**_Accept
*****0.01_****g?_**_Under
    
```

A.Print: off; Stable: off  
 Stable reading,  
 Unstable reading

**Scout Pro Print Format 1 (for the models 303/123/202/402/602/2001/6001/401FZH/601FZH/6000FZH):**

Output String:

[weight]	12 characters (right justified)
[space]	1 character
[unit]	5 characters (left justified)
[space]	1 character
[stability indicator]	1 character; "?" when unstable, space when stable
[Legend]	1~10 characters

Output example:

```

*****0.00_g****_*
*****12.73_g****_?
*****0.85_oz***_WET*WT

```

**Scout Pro Print Format 2 (for the models 401/601/6000):**

Output String:

[weight]	11 or 12 characters (right justified)
[space]	1 character
[unit]	1~5 characters
[space]	1 character
[stability indicator]	1 character; "?" when unstable, space when stable
[space]	1 character
[Legend]	1~10 characters

**Note:** The unit field length varies with different units. The weight field could be 11 or 12, depending on if the weight string has a dot or not.

Output example:

```

*****100_g_*
*****273_g_?_
*****8.5_oz***_WET*WT

```

**Print Format 3 (for Certain POS Systems):**Output String:

[weight]	11 characters (right justified)
[space]	1 character
[unit]	5 characters (right justified)
[stability indicator]	1 character; "?" when unstable, space when stable
[Term]	2 characters

**Note:** All of the fields have fixed length.

Output example:

```
*****0.00_****g*  
*****12.73_****g?
```

## USB DEVICE INPUT

The following interface commands will be acknowledged by the Balance. They are case sensitive. The balance will return “ES” for invalid commands.

Command	Function
IP	Immediate Print of displayed weight (stable or unstable).
P	Print displayed weight (stable or unstable).
CP	Continuous Print.
SP	Print on Stability.
SLP	Auto Print stable non-zero displayed weight.
SLZP	Auto Print stable non-zero weight and stable zero reading.
xP	Interval Print x = Print Interval (1-3600 sec) 0P ends interval Print
0P	Turn off Auto Print
H	<b>H x “text”</b> Enter Header line , where x = line number 1 to 5, “text” = header text up to 24 alphanumeric characters
F	<b>F x “text”</b> Enter Footer line , where x = line number 1 to 2, “text” = footer text up to 24 alphanumeric characters
Z	Same as pressing Zero Key
T	Same as pressing Tare Key.
xT	Establish a preset Tare value in displayed unit. x = preset tare value. Sending 0T clears tare (if allowed).
PT	Prints Tare weight stored in memory.
PM	Print current application mode (weighing mode).
xM	Set current application mode to x. x depends on application 1M: WEIGH, 2M: COUNT, 3M: PERCENT, 4M: CHECK, 5M: DYNAMIC, 6M: TOTAL, 7M: DENSITY, 8M: HOLD, 9M: MOLE
M	Scroll to the next enabled mode.
PU	Print Current weighing unit: g, kg, lb, oz, etc....
xU	Set balance to unit x: g, kg etc. 1U: g, 2U: kg, 3U: ct, 4U: N, 5U: oz, 6U: ozt, 7U: dwt, 8U: lb, 9U: lb:oz, 10U: gm, 11U: thk, 12U: tsg, 13U: ttw , 14U: tola, 15U: c
U	Scroll to the next enabled unit.
ON	Brings out of Standby
OFF	Goes to Standby.
C	Begin Span Calibration
AC	Abort Calibration.
PSN	Print Serial Number.
PV	Print Version: print name, software revision and LFT ON (if LFT is set ON).
x#	Set Counting APW (x) in grams. (must have APW stored)
P#	Print Counting application APW.
x%	Set Percent application reference weight (x) in grams. (must have reference weight stored)
P%	Print Percent application reference weight.
xCO	Set Checkweighing Over Limit in grams x.
xCU	Set Checkweighing Under Limit in grams x.
PCO	Print Checkweighing Over Limit.
PCU	Print Checkweighing Under Limit.
xMM	Set Molar Mass in g/mol .
PMM	Print Molar Mass
xS	0 = print unstable data, 1 = print stable only
xFMT	0 = New Scout print format (default) ; 1 = Scout Pro print format 1; 2 = Scout Pro print format 2; 3 = for certain POS system.
xRL	0 = disable response; 1 = enable response.

## ACCESSORIES

For a complete listing of OHAUS printers and other accessories, contact OHAUS Corporation or visit [www.ohaus.com](http://www.ohaus.com).

## COMPLIANCE

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Please note that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## DISPOSAL

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related. For disposal instructions in Europe, refer to [www.OHAUS.com/weee](http://www.OHAUS.com/weee). Thank you for your contribution to environmental protection.

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



Menu supplémentaire	Fonction
Communication  --Clef USB  --Débit en bauds	Réglez la vitesse de transmission (bits par seconde).  <p>The screenshot shows a 'Baud Rate' menu with a grid of options: 1200, 2400, 4800, 9600, 19200, 38400, 57600, and 115200. The '9600' option is highlighted. Below the grid are 'Back' and 'Exit' buttons.</p>
--Transmission	Réglez les bits de données, le bit d'arrêt et la parité.  <p>The screenshot shows a 'Transmission' menu with a grid of options: 7E1, 7E2, 7N1, 7N2, 7O1, 7O2, 8N1, and 8N2. The '8N1' option is highlighted. Below the grid are 'Back' and 'Exit' buttons.</p>
--Établissement de liaison	Définissez la méthode de contrôle de flux.  <p>The screenshot shows an 'RS-232' menu with three options: 'Baud Rate' (with a clock icon), 'Transmission' (with a data icon), and 'Handshake' (with a handshake icon). The 'Handshake' option is highlighted. Below the grid are 'None' and 'Xon/Xoff' buttons.</p>
--Paramètres de l'impression	Pour davantage d'information, se reporter au mode d'emploi STX.  <p>The screenshot shows a 'Communication' menu with two options: 'RS-232' (with a printer icon) and 'Print Settings' (with a printer icon). Below the grid are 'Back' and 'Exit' buttons.</p>

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<p><b>PRINT</b>                  Reset: <b>no</b>, yes                  Stable: <b>off</b>, on                  A.Print:                      <b>on</b>.Stable                      interval                      accept                      continuous                  Content: Result (-&gt; <b>off</b>, <b>on</b>)                      Gross (-&gt; <b>off</b>, <b>on</b>)                      Net (-&gt; <b>off</b>, <b>on</b>)                      Tare (-&gt; <b>off</b>, <b>on</b>)                      Header (-&gt; <b>off</b>, <b>on</b>)                      Footer (-&gt; <b>off</b>, <b>on</b>)                      Mode (-&gt; <b>off</b>, <b>on</b>)                      Unit (-&gt; <b>off</b>, <b>on</b>)                      Info (-&gt; <b>off</b>, <b>on</b>)                  Layout: Format (-&gt;Single,<b>Multi</b>)                      Feed (-&gt;<b>Line</b>, 4 Lines, Form)                  End Print: Exit menu</p>	<p><b>USB.Dev</b>                  Reset: <b>no</b>, yes                  Baud Rate: 1200...<b>9600</b>...115200                  Parity: 7 even, 7 odd, 7 none,                      <b>8 none</b>                  Stop bit: 1, 2                  Handshake: <b>none</b>, XOn-XOff                  Alternate command: Print (A...P...Z)                      Tare (A...T...Z)                      Zero (A...Z)                  End USB.Dev: Exit menu</p>
--	---

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P/N 30269075 C © 2016 Ohaus Corporation, all rights reserved / todos los derechos reservados/ tous droits réservés/ Alle Rechte vorbehalten / Tutti i diritti riservati / 版权所有 / 모든 권리 보유 / Все права защищены / すべての権利予約