

VETEK



**Personal scale WBA-300 with
PROFESSIONAL ELECTRONIC COMPOSITION**

USER MANUAL



Read this manual carefully before using the instrument

INDEX

2.	GENERAL RULES.....	PAG. 2
3.	SAFETY.....	PAG. 3
4.	THECHNICAL SPECIFICATION.....	PAG. 4
5.	INDICATOR.....	PAG. 5
6.	INSTALLATION.....	PAG. 7
7.	GENERAL CONDITIONS FOR CORRECT MEASUREMENTE.....	PAG. 8
8.	INSTRUCTION FOR USE.....	PAG. 8
9.	PRINT FUNCTION.....	PAG. 12
10.	MAINTENANCE AND ASSISTANCE.....	PAG. 13
11.	COMFORMITY.....	PAG. 13
12.	SCRAPPING AND WASTE DISPOSAL.....	PAG. 13
13.	WARRANTY.....	PAG. 14
14.	THERMAL PRINTER.....	PAG. 15
15.	ASSEMBLY-USE MECHANICAL TELESCOPIC HEIGHT METER (OPTIONAL).....	PAG. 18

By choosing the **WUNDER** mod. **WBA** professional electronic scale, you have purchased a high precision instrument.

Since over 40 years **WUNDER** has placed its experience at the service of health. This instrument is calibrated in conformity with accuracy class III.

The instrument is characterized by the possibility of fixing the electronic terminal to the weighing platform by means of a column or else installing the terminal autonomously.

The instrument is equipped with a dual LCD electronic terminal with triple reading to simultaneously view Weight, Height, and BMI

2. GENERAL RULES

 **Carefully read this manual before using the instrument as it supplies**

important indications concerning OPERATING SAFETY AND MAINTENANCE.

The descriptions and illustrations supplied in this manual are not binding.

WUNDER reserves the right to modify its products as deemed convenient in order to approve them without being committed to update this publication.

Conventions:

The following symbols have been used in this manual:



0470 MEDICAL CERTIFICATION



METROLOGICAL CERTIFICATION AND APPROVAL FOR MEDICAL USE



ACCURACY CLASS



ASSISTANCE WARRANTY CERTIFICATE



ATTENTION!

PLACED BEFORE DETERMINING PROCEDURES. COMPLIANCE FAILURE CAN DAMAGE THE PRODUCT

WARNING!



PLACED BEFORE DETERMINING PROCEDURES. COMPLIANCE FAILURE CAN HARM THE OPERATOR OR PATIENT OR DAMAGE THE PRODUCT



WASTE DISPOSAL



TYPE B PARTS SUPPLIED



BATTERY POWER



INDICATION OF WEIGHT FUNCTIONALITY

INDICATION OF STABLE WEIGHT



INTERFERENCE COULD OCCUR NEAR THE APPLIANCES



DUAL INSULATION

3. SAFETY



WARNINGS!

Operators must read this manual carefully, comply with the instructions it contains and become familiar with the correct use and maintenance procedures of the instrument.

- Retain this manual for consultation and as a help in staff training
- Do not overload the instrument beyond its maximum capacity
- Do not apply loads abruptly.
- Do not press the keys with sharp or pointed objects
- Do not try to open the instrument.
- Do not remove seals from the instrument.
- Do not short-circuit the battery terminals
- Use only the power supply provided by Wunder. Before using it, make sure that the local mains voltage is compatible with the voltage of the adapter shown on the identification plate
- Regularly check the integrity of the instrument's power cord and make sure it does not come in contact with hot appliances
- Make sure that the power cord does not create obstruction hazards
- Unplug the instrument before cleaning it
- Do not place the instrument in water or other liquids
- Perform maintenance and subsequent metric verifications regularly



WARNING : People with pacemakers and other medical devices internal.

Note that people with an electronic medical implant, such as a pacemaker, should not use WBA, as it passes a low-level electrical signal through the body, which may interfere with its operation.



WARNING!

For the use of the instrument by persons with disabilities, it is recommended to include qualified assistance.

Keep the instrument out of the reach of children.



WARNING:

The body composition analyzer, should be used with bare feet. Be sure to clean the weighing platform with appropriate disinfectant after each use.

Do not pour any liquid directly on the weighing platform, to avoid damaging it.

To clean the platform, use a soft cloth moistened with ethyl alcohol.


Avoid the use of harsh chemicals, excessive amounts of water or high pressure washing systems.

When cleaning, always remove the power cable, avoid touching it with wet hands.

4. THECHNICAL SPECIFICATION

Model	WBA	
Capacity	300g	
Division	100g	
Unit weight	kg	
OIML	Class III	
Display	1" LCD weight 5 digits 20mm ; 2" LCD BMI/ % Fat Mass	
Function Keys	ON-OFF, Zero, Tare, Printer; Body type; Sex, Time, Numerical keys 0~9; Reset key ; Enter key	
Technology BIA	Bioelectrical Impedance Analysis 4 feet- electrode	
Impedance frequency	50kHz 500uA	
Impedance Range	200 ~ 1000Ω / 0.1Ω	
Age Range	10 ~ 80 years, increase of 1 year	
Height Range	60 ~ 210cm / 3ft ~ 7ft 11.0in	
Power	Rechargeable Batteries 7.2V 2000mA o batteries AA (1.5X6) 9V.Medical Adapter AC 12V 2A	
Platform dimension	340 (L) x 450 (P) x 80 (H) mm	
Temperature Use	5°C ~ 35°C - 30% / 80% UR	
Storage temperature	-20°C ~ 60°C - 0% / 90 % UR	
Output	RS232 9 pin female	
Optional (not included)	External thermal printer; Telescopic height meter	
Printer Data	Date	GG/MM/AAAA hh:mmEs: 30/10/2013 10:55
	Mode	Standard / Athletic
	Sex	Male / Female
	Age	10 ~ 80 years, increase of 1 year
	Height	60 ~ 210cm / 3ft ~ 7ft 11.0in
	Weight	0 ~ 300 kg
	BMI	Increase of 0.1
	BMR	Increase of 1 kcal
	Fat %	5 ~ 50% / increase of 0.1%
	Fat kg	Increase of 0.1 kg
	FFM	Increase of 0.1 kg
	TBW	Increase of 0.1 kg

5. INDICATOR



1. Specifications
2. Main functions
3. Numeric Keyboard
4. Display Weight
5. Display BMI / % Fat Mass
6. Power Indicator

Keyboard

Key	Description
	On/Off key. Pressed for 3 seconds it switches off the scale.
	Indication reset (within ±2% of the capacity)
	Printer Data
	Sex female
	Sex Male
	ATHLETIC MODE: "means a person subject to intense physical activity for at least 12 hours a week and with a resting heart rate of 60 beats per minute or less. They are also considered athletic people if for some years they have practiced sport and they currently exercise less than 12 hours per week. Are not considered in this mode, beginners while exercising for at least 12 hours a week have not yet taken a defined athletic body shape.
	Time setup
	Tare Function of unwanted weight (clothes, etc ...). Function net weight.
	Data Cancel
	ENTER key
	Numeric keyboard for entering data

Display functions:

DATE Mark
To set date and time.

Power Index
100% → 75% → 50% → 25%

Athlete Mark
Indicates that the user is an "athlete"

BMI
Body Mass Index is a height to weight ratio, and is calculated by the following formula

Body fat percentage
Indicates that the figure shown is the body fat percentage

PRE- TARE HOLD NET WEIGHT

lb

kg

ATHLETIC

BMI

FAT(%)

cm

AGE

STEP ON

PRE- TARE HOLD NET WEIGHT

Unit Mark (kg)
Indicator that the weight is shown in kilograms

Gender Mark (♀ / ♂)
The user is male or female.

Unit of height (cm)
Indicates that the height is shown in centimeters

Age Mark
User's age

Ready Mark
The device is ready to begin measurements.

kg

cm

AGE

STEP ON

ATHLETIC

BMI

FAT(%)

6. INSTALLATION



AFTER HAVING VERIFIED THE INTEGRITY OF THE EQUIPMENT AND OF ALL ITS COMPONENTS, CONNECT THE INSTRUMENT WITH THE

SUPPLIED POWER SUPPLY TO CORRECTLY RECHARGE THE BATTERIES

Position the instrument on a level and stable surface

The instrument must be placed away from heat sources, in an environment without excessive vibrations and drafts.


Verification of equipment:

Model WBA 300 with printer

1. Base scale with support for column
2. Column and display
3. Thermal printer mounted on the column with a metal support with of the printer: already mounted on the column
4. External power supply for the scale
5. External power supply for the printer
6. User manual
7. Height meter (optional not included)

Assembly (for column version only)

Figure A



Assembly is easy: fix scale's column to the weighing platform (See Figure A)



ATTENTION!

FOLLOW THE INDICATIONS BELOW BEFORE INSERTING THE COLUMN

Figure B





Figure C



- Overturn the weighing platform resting it onto a level surface.
- Remove the display connector from inside the column and connect it to the weighing platform connector, making it pass through the hole indicated by the red arrow (see Figure B).
- Fix the column to the weighing platform support using the 4 supplied screws. (see Figure C).
- Position the instrument on the floor and adjust the 5 feet to it.
- To use, connect the supplied external power supply to the current (see battery operation)

7. GENERAL CONDITIONS FOR CORRECT MEASUREMENT



WARNING!

Preventive steps to be carried out in advance of the test execution

1. To ensure accuracy, readings should be taken away from the meals. This is necessary since the mass of the food is measured together with the weight and it is susceptible to precise measurement.
2. Prior to testing, if possible, you should urinate and defecate. Although not included in the elements that constitute the body mass, the volume of urine and excrement is included in the measurement.
3. Did not run in conjunction with the following medical devices:
 - Pacemakers and electro medical equipment like (Bypass).
 - Systems of artificial life, such as artificial hearts or lungs
 - Portable electrical devices, such as electro-cardiographs
4. Trying to be at rest before performing the test: strenuous exercise or sudden movements can cause temporary changes in body composition.
5. Whenever possible, avoid performing the test after taking a shower or a sauna
The exudation causes temporary changes in body composition.
6. Did not take the test during the menstrual cycle
Women experience an increase in water in the body during the menstrual cycle
7. Perform the test at normal temperatures (20 ° C-25 ° C)
While the human body is stable at normal temperatures, body composition is susceptible to extreme heat or cold.
8. Always remember that more the body condition of the various measurements are similar and standards (eg wearing the same clothing, running the test before eating and do any type of exercise, etc. ..), plus the test results are reproducible and reliable.


8. INSTRUCTION FOR USE

8.1 Instruction for charging and connecting

If  prompt displays on the LCD, please charge the scale with exclusive adaptor.



Locate adaptor plug on the rear side of indicator.

The battery should be recharged at least every 3 months regardless if it is used or not. After a long period in storage, e.g. over 3 months, the battery should run a full cycle (charge/discharge) to allow it to restore to full capacity.



Ac Jack connector

8.1a Rechargeable Batteries







1. Open the battery cover, and take the battery out.
2. Disconnect the rechargeable battery pack then replace it by a new one.
Put the battery back then arrange the connector in place well.

8.2 DATE AND TIME SETTING

Keep the HOLD key pressed for 3 seconds to access TIME SETTING programming.

Example: to program 25 December 2008, 8:00 a.m.

	Year programming: use the number keys to enter the correct value, press HOLD to access the next step once correctly set.
	Date programming: Use the number keys to enter the correct value, press HOLD to access the next step once correctly set.
	Time programming: Use the number keys to enter the correct value, press HOLD to confirm the entered date and time once correctly set.
	Display format: YYYY → MM.DD → HH:MM

8.3. ONLY WEIGHT:

Questo strumento può essere utilizzato anche come normale bilancia.


1. Press the power button 
2. Place the patient and the first display will show the weight.


8.4 WEIGHT & BODY COMPOSITION

 **WARNING!**


THE INFORMATION PROVIDED BY THE SCALE AND USED TO DETERMINE DIETS AND EXERCISE PROGRAMS,
MUST BE READ AND CONSTRUED BY A QUALIFIED





1




Press the key  the top display shows 0.0kg


2




Select the body conformation with the key  or,  OR  and press the key  to confirm.


3




Insert height and press 

4



Insert age eight and press 

5




A flashing arrow appears on STEP ON.

Walk up barefoot on the 4-electrode of the platform and stand still

The upper display shows the weight detected (Figure 5).


5



The lower display shows 'oooo' and begins the analysis of the body.


When the symbols 'oooo' is moving, indicate the measurement body of the analysis and...

6



.....the lower display shows BMI and %

7



Press  to displays the age and height



To print the results control that the printer is on

Press the key  of the printer to print the results.

The measured data can be printed using thermal printer: date and time; Data included: body type; Data measured: weight, BMI, Body Fat%, Fat Mass kg, kg lean body mass, body fluids, Basal Metabolism and impedance.


8.5 TARE (CLOTHES):

This feature allows you to remove the weight of the clothes from the measurement of the weight. Turn the instrument on by pressing .

When the display shows  and you hear a beep, press the key . The display shows 0000 flashing.

Select with the numeric keyboard the value of the cloths **ES. 001,5kg**. To memorize the weight press the key .

The display will show the selected value in the negative. The value remains in memory until it is changed or deleted.


To delete the value PRE-TARE repeat the process by selecting the numeric keypad and pressing 0000 then press  confirm.

9. PRINT FUNCTION


Standard	EIA RS-232 Compatible
Method of communication	Transmission asynchrony
Speed	9600 bps
Data length	8 bits
Parity Check	None
Stop bit	1 bit
Handshake	None
Code	ASCII

Linea serial RS232

Terminal	Signal
2	TX
3	RX
5	GND



In case of conection to PC or a printer TP2100, use the serial cable invert (pin 2 e 3).



PRINT RESULTS

WUNDER	
15/ 03 / 2013 10: 30	
Weight	70,5kg
Height	180,0cm
Tare: weight clothes	1,5kg
BMI: Body Mass Index	25,3
FAT: Fat Mass	25,5%
FAT: Fat Mass	15,7kg
FMM: Free Mass	55,4kg
BMR: Basal Metabolic	1660 ckal
TBW: Total Body Water	41,8kg
Impedance	603Ω

The format presented below is the standard format of results print-out and cannot be changed.

More information info@wunder.it or service@wunder.it

Start Menu → Programs → Accessories → Communication → Hyper Terminal.

- New Connection Description

Give new connection a name then click OK.

- Select Your COM Port

Click Connect to select your COM port. Usually there's only one option for select. Then click OK.

- Port Settings

Click Bits per second to set up rate at 9600, Data bits at 8, Parity at None, Stop bits at 1 and Flow control at Hardware. Then click OK to complete your setting.

- Output Data

When press the PRINT key to output data from scale to PC or an Optional Printer.

10. MAINTENANCE AND ASSISTANCE

For better and longer duration of the product it should receive thorough general cleaning periodically.

The instrument must be cleaned with a soft cloth moistened with water or neutral detergent, without using solvents or abrasive substances. If the instrument remains idle for a long period, remove the batteries from the terminal.

During shipping, make sure not to subject the instrument to blows or excessive mechanical stress.

In case of repairs or assistance, contact your dealer or an authorized service centre contacting service2@wunder.it or sales@wunder.it.

11. CONFORMITY

WUNDER MODEL WBA CABLE ELECTRONIC BODT FAT COMPOSITION SCALE SERIAL N°

We hereby certify that this instrument has been inspected and has successfully passed the functional test.
It complies with the following standards and directives:

**EN 45501 / EN60601-1-2 / EN60601-1
93/42/EEC – Medical Devices Directive**

12. SCRAPPING AND WASTE DISPOSAL

If set aside for a long period, protect those parts which could be damaged due to dust build-up

Scrapping

When you decide to no longer use this item, we recommend making it unusable. We also recommend making those parts which could be sources of danger harmless

Waste disposal  EU 2002/96/EC

This product complies with the EU Directive 2002/96/EC. The symbol of the crossed-out waste bin on the appliance indicates that the product, needing to be treated separately from household waste, at the end of its useful life must be completed in a separate collection facility for electric and electronic appliances or returned to the dealer upon purchase of a new equivalent appliance. The user is responsible for bringing the appliance to an appropriate collection structure at the end of its life. Appropriate separate collection and sending the appliance for recycling, treatment and environmentally compatible waste disposal contributes to avoid possible negative effects on the environment and health and favours the recycling of the materials the product is made of.

For more detailed information regarding available collection systems, contact your local waste disposal service or the shop where the product was purchased.

As consumers, you are obliged by law to return used or dead batteries. You may deposit old batteries at public collection spots in your town or else with any battery dealer who has placed specific collectors for this purpose. Even when scrapping electric and electronic appliances, they must be removed and deposited in specific collectors.

NOTE: The following symbols indicate the presence of harmful substances

Pb Pb = batteries containing Lead

Cd Cd = batteries containing Cadmium

Hg Hg = Batteries containing Mercury



ATTENTION!

Do not throw electric parts and used batteries away with household waste. Dispose of the batteries by means of your closest collection centres.

13. WARRANTY

This certificate must be kept until the warranty has expired.

It must be presented together with the invoice, tax receipt or delivery note providing the name of the dealer and date of purchase whenever a technical intervention is required. Otherwise the user will lose any warranty rights.

The warranty takes effect from the date of purchase and is valid during the entire period foreseen by the current catalogue/pricelist.

By warranty we mean the replacement or repair free of charge of parts making up the appliance which, at the discretion of the manufacturer, are deemed faulty from the origin; Wunder therefore has the faculty of repairing or replacing the item.

The warranty does not cover:

- shipping faults, damage caused by falls, carelessness or tampering
- damage caused by incapability of using the appliance and of its improper use
- damage caused by an insufficient or inadequate electrical system or alterations resulting from environmental, climatic or other types of conditions
- damage due to incorrect installation of the appliance and repairs carried out by unauthorized personnel
- Interventions at home for convenience controls or presumed defects
- Routine maintenance and that which can be considered normal wear from use
- consumables such as: power supplies, batteries, keyboards, plates, wheels, heads, rolls, load cells faulty due to blows or overloads

Service can also be refused when the appliance has been changed or transformed in any way.

In case of interventions at one's home, the customer must pay the fixed fee; if however the appliance is repaired at an authorized Wunder Service Centre, expenses and relative travel risk are borne by the user.


Wunder will not be held liable for damage of any nature caused directly or indirectly to persons, animals or objects resulting from failure to comply with all the instructions indicated in this manual or anyway resulting from improper use.

The Court of Bergamo has jurisdiction in case of any dispute

14. THERMAL PRINTER (OPTIONAL)




 **Read this manual carefully before using the instrument**



14.1. HOW TO USE YOUR PRINTER


Connect the AC Adaptor and RS232 interface cable. Press the I/O button to power on the printer.



The PRINTER XL printer can be connected via the RS232 port and cable with the RA200 XL.


Press the PRINT key.

Print double language


**CHANGING THE PAPER ROLL**

The paper print roll on a PRINTER XL printer is changed as follows:

1. Open the upper cover and position the paper roll so that it feeds in the right direction.



2. Feed the end of the paper roll through the printer cover and close the cover. The printer is ready to print.



--	--

14.2. PRINTER SPECIFICATION

ITEM	SPECIFICATION
Printing method	Thermal dot line printing
Number of dots	384 dot/line
Resolution (dots/mm)	8 dot/mm
Printing width (mm)	48 mm
Max. paper thickness (μm)	80 μm
Paper	57.5 \pm 0.5 mm (W) x ψ 30mm max
Paper loading	Easy loading
Printing speed (mm/s)	75 mm/s
Dimensions (L x W x H)	120 x 80 x 50mm
Operating temperature ($^{\circ}\text{C}$)	0 $^{\circ}$ to 50 $^{\circ}\text{C}$
Operating humidity (%)	20% to 85% (no condensation)
Interface	Conform to RS232
Voltage range (adaptor only)	12V 2A

SHOCK HAZARD WARNING

The printer and power supply should never be operated in a location where either one can get wet. Personal injury could result.

MEDIA WARNING

Always use high quality approved paper print roll. The approved product can be ordered from your dealer.

RELOADING HINT

If printer runs out of paper roll while printing, switch the power off and reload paper roll. Power on and print again.

14.3. ERROR MESSAGES

Problem	Solution or Reason
Status Indicator does not light when power switch is turned to ON (I) position.	Check the power connections from AC outlet to the printer power supply.
Status indicator illuminated (green light), but the printer will not print.	Check interface cable connections from the scale to printer. Make sure the top cover is closed and locked.
Printer paper jammed.	Reset the paper roll. Please refer to "CHANGING THE PAPER ROLL".
Status Indicator is flashing green.	Paper roll is out. Reload a new source of paper roll so the printer can continue printing. Please refer to "CHANGING THE PAPER ROLL".

15. ASSEMBLY-USE MECHANICAL TELESCOPIC HEIGHT METER (OPTIONAL)










ATTENTION!


ONLY FOR COLUMN SCALE

HEIGHT ROD KIT


1. Metal support for telescopic height meter		 side view metal support
2. Black abs support height meter		

3.Screws	 A)	 B)	 C)	 D)	 E)
	N°02 screws T.E. (M5x16)	N°02 slotted screws (M5x25)	N°01 Allen Screws (M6x12)	N°02 nuts \varnothing 5	N° 02 washers \varnothing 5
4.Tools					
	Tools		Allen tools		


Picture 1: Assemble white metal support with abs black part with the screws (B+E). Lock all with the screws (D)










Picture 2-2a: behind the pole (column) assemble the white metal support with screws (B). **Picture 2b:** the position is on the right



Picture 3-3a: assemble the height rod on black abs support with screws (C) and allen tools kit **Picture 3b:** The height rod is ready for use.



MEASURING OVER 131 cm		MEASURING UNDER 130 cm	
	Pull up the height meter rod		Press the black button and measuring height and push down

 <p>2</p>	<p>Position it horizontally</p>	 <p>2</p>	
 <p>3</p>	<p>Lower the rod up to the point where the lever rests head touches the patient's head. The height can be read at the point 'A'</p>	 <p>3</p>	<p>Position it horizontally</p>
		 <p>4</p>	<p>Lower the rod up to the point in which the blade rests head touches the patient's head. The height can be read at point 'B'.</p>