



www.wunder.it

PROFESSIONAL ELECTRONIC BABY SCALE

MOD. BABY02

USER MANUAL





Read this manual carefully before using the instrument

2. GENERAL RULES	2
3. SAFETY	3
4. SPECIFICATION	6
5. INDICATOR	6
6. INSTALLATION	7
7. ISTRUCTION FOR USE	8
8. PRINTER FUNCTIONS	11
9. ERROR MESSAGE	12
10. MAINTENANCE AND ASSISTANCE	
11. CONFORMITY	13
12. IDENTIFICATION LABELS	
13. SCRAPPING AND WASTE DISPOSAL	13
14. WARRANTY	14
15. INSTALLATION BABY HEIGHT ROD	15
16. INSTALL SECURITY BRACKETS	16

By choosing the WUNDER mod. Baby02 professional electronic scale, you have purchased a high precision instrument.

Since over 40 years **Wunder** has placed its experience at the service of health. This instrument is compliant with national standards in hospitals and clinics with medical class I with measurement function and is calibrated in conformity with accuracy class III.

The instrument is characterised by the possibility of fixing the electronic terminal to the weighing platform by means of a column or else installing the terminal autonomously.

The instrument is equipped with a dual LCD electronic terminal with triple reading to simultaneously view Weight and weight difference calculation.

2. GENERAL RULES



Carefully read this manual before using the instrument as it supplies

important indications concerning OPERATING SAFETY AND MAINTENANCE.

The descriptions and illustrations supplied in this manual are not binding.

WUNDER reserves the right to modify its products as deemed convenient in order to approve them without being committed to update this publication.

Conventions: The following symbols have been used in this manual:



0476 EC TYPE APROVAL FOR MEDICAL USE



METROLOGICAL CERTIFICATION AND APPROVAL



ACCURACY CLASS





PLACED BEFORE DETERMINING PROCEDURES. COMPLIANCE FAILURE CAN HARM THE OPERATOR OR PATIENT OR DAMAGE THE PRODUCT



WASTE DISPOSAL



TYPE B PARTS SUPPLIED



BATTERY POWER



INDICATION OF WEIGHT FUNCTIONALITY



INDICATION OF STABLE WEIGHT



INTERFERENCE COULD OCCUR NEAR THE APPLIANCES



DUAL INSULATION (CLASS II)

3. SAFETY



ATTENTION!

Operators must read this manual carefully, comply with the instructions it contains and become familiar with the correct use and maintenance procedures of the instrument.

- Retain this manual for consultation and as a help in staff training
- Do not overload the instrument beyond its maximum capacity
- Do not apply loads abruptly.
- Do not press the keys with sharp or pointed objects
- Do not try to open the instrument.
- Do not remove seals from the instrument.
- Do not short-circuit the battery terminals
- Use only the power supply provided by Wunder. Before using it, make sure that the local mains voltage is compatible with the voltage of the adapter shown on the identification plate
- Regularly check the integrity of the instrument's power cord and make sure it does not come in contact with hot appliances
- Make sure that the power cord does not create obstruction hazards
- Unplug the instrument before cleaning it
- Do not place the instrument in water or other liquids
- Perform maintenance and subsequent metric verifications regularly

3.1 Manufacturer's guide and declaration – electromagnetic immunity

Guide and Statement of manufacturer - Electromagnetic emissions

The electronic scales **BABY 02** model is scheduled for operation in the electromagnetic environment specified below.

the customer and the user should ensure that it is used in that environment

Emission test	Conformity	Guide to electromagnetic environment	
		The baby scales BABY 02 model uses RF energy only for its	
RF Emission CISPR11	Group 1 Class B	internal function. The refore its RF emissions are very low	
NI EIIII331011 CI31 NIII		and probably not cause no interference in electrical	
		equipment.	
RF Emission CISPR11	Group 1, Class B	The heliconder BADY 03 is suitable for one in all conitations	
Harmonic emission	Class A	The baby scales BABY 02 is suitable for use in all sanitary	
Voltage fluctuations/		and hospital buildings, connected to the low voltage public	
flicker emission	Compliant	power supply network.	

Guide and Statement of manufacturer-Electromagnetic Immunity

The electronic scales BABY 02 model is scheduled for operation in the electromagnetic environment specified below.

the customer and the user should ensure that it is used in that environment.

Immunity test	Compliance	Guide to electromagnetic environment
Electrostatic discharges (ESD) IEC/EN61000 - 4 -2	6kV in contact 8kV in air	The floors should be made of wood, concrete or ceramic. If the floors are covered in synthetic material, the relative humidity should be maximum 30%.
Rapid transients/burst IEC/EN61000 - 4 - 4	2kV power supply	The power supply should be of the type used typically in commercial or hospital environments.
Surge IEC/EN61000 - 4 - 5	1kV differential mode	The power supply should be of the type used typically in commercial or hospital environments.
Voltage dips, short interruption and voltage variation IEC/EN61000 - 4 - 11	<5%UT for 0.5 cycle 40%UT for 05 cycle 70%UT for 25 cycle <5%UT for 5 sec	The power supply should be of the type used typically in commercial or hospital environments. Note=UT is the value of the voltage of the feed.
Power frequency Magnetic field IEC/EN61000 - 4 - 8	3A/m	

NOTE: UT is the value of the power supply voltage.



CAUTION

The medical device requires particular electromagnetic compatibility precautions and must be installed and used according to the information provided in the accompanying documents.

Manufacturer's guide and declaration - Electromagnetic emissions			
The electronic scales BABY 02 model can be used in the electromagnetic environment specified below. The			
customer or the user	customer or the user of the seat must ensure that it is used in compliance with the specifications below.		
Emissions test Compliance Guide to electromagnetic environment		Guide to electromagnetic environment	
Radiated/	3Vrms 150kHz to	Equipment in RF communications and mobile laptops should ot	
conducted / 80MHz(for be used near any part of the tooth, including cables, except		be used near any part of the tooth, including cables, except	
emissions appliances that are when they meet the recommended separation distance		when they meet the recommended separation distance	
IEC/EN61000-4-6 not life supporting) calculated by the frequency applied to the transmitter.			

		Recommended separation distance d = 1.2 V P d = 1.2 V P from
		80 MHz to 800 MHz d = 2.3 V P from 800 MHz to 2.5 GHz where
		the P and maximum rated power output of the transmitter in
 Radiated/	3Vrms 150kHz to	Watts (W) according to the manufacturer of the transmitter and
		the recommended separation distance in meters (m).
conducted	80MHz(for	The intensity of the field of a fixed RF transmitters, as
/emissions	appliances that are	determined by an electromagnetic site, it could be lower than
IEC/EN61000-4-3	not life equipment)	the level of compliance of each range of frequency b. You can
		check interference in proximity of equipment marked with the
		following symbol:

With 80 MHz and 800 MHz is applied the higher frequency range. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) The intensity of the field for fixed transmitters such as base stations for radiotelephone mobile and cordless phones) and land radio mobile, amateur radio equipment, radio transmitters in the AM and FM and TV transmitters cannot be made theoretically and with precision. To establish an electromagnetic environment due to fixed RF transmitters, it should consider the electromagnetic survey of the site. If the field strength measured at the place where you use the instrument exceeds the applicable level of compliance of the above, it should be placed under observation normal operation of the tooth divice. If you notice abnormal performance, it may take additional measures as a different orientation or position of instrument.

of frequencies from 150 kHz to 80 MHz should be less than 3 V / m.

Recommended separation distances between portable radio equipment unity.

The scales **BABY 02** model is expected to operate in an electromagnetic environment in which they control the noise radiated RF. The customer or the tooth may be of help to prevent electromagnetic interference by ensuring a minimum distance between mobile communication devices and portable RF (transmitters) and unity, as recommended below, with respect to the maximum power output of equipment radio.

Output power rating	Distance separating the frequency of the transmitter m		
of the transmitter W	150 kHz a 80 MHz d =	80 MHz a 800 MHz d =	800 MHz a 2,5 GHz d =
of the transmitter vv	1,2 √P	1,2 √P	2,3 √P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with maximum rated power output not reported above, the recommended separation distance d in meters (m) can be calculated using the equation applicable to the frequency of the transmitter, where P and the maximum rated power output of the transmitter in Watt (W) according to the manufacturer of the transmitter.


Note: With 80 MHz and 800 MHz is applied the higher frequency range.

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

4. SPECIFICATION

Model	BABY 02	
Capacity	Max 3/6kg, e=1/2g; Max 6/15kg, e=2/5g; Max 10/20kg, e=5/10g;	
Approvazione OIML	Class III	
Unit weight	kg	
Display	Double LCD: 1 LCD 20mm weight with 5 active digits; 2 LCD 12mm	
Dimensions	550 x 290 (DISH)	
Key functions	ON/ZERO/OFF, HOLD, TARA, PRINT, MEMORY(MILK BOTTLE)	
Power	Adapter 12 V - 2 A Rechargeable Batteries Pack- (OPTIONAL)	
Temperature operation	5°C / 35°C	

5. INDICATOR



- 1. Weighing capacity and division
- 2.Key function
- 3.Main LCD: Weight
- ${\it 4.Second\ LCD\ automatically\ calculates\ the}$
- amount of Milk consumed
- 5. Power display

Key functions

Tasto (b) →O←

Description

Key ON. To reset the display (OFF) long press 3 seconds. Zero the scale (±2% of full capacity)



To hold weighing result on the display/ determine stable weighing value.



Tare un-wanted weight.




Function Weight-Milk: automatically calculates the amount of Milk consumed



Print out data.

Function Display Symbols



Stable symbol: To indicate that the weight is stable.

Minus weight: Weight under zero.

Zero symbols: Weight is at zero point.

Low battery: Battery need to charge or replaced.

6. INSTALLATION

After removing the instruments from the packing check the integrity and right contents Put the scale far from source from of heat e vibration.

Put the scale on stable surface.

Indicator assembly (is not assembled)

Run the cable through the tube





Run the cable through the tube. Setup the bracket. Screw to the holes.



Indicator Column assembly with platform








1. Run the cable through the column.


2. Connect the jack

3. Tighten the three screws of the base by turning

cable



4. Align the platform on the base and tighten the screws



5. Fit the adjustable foot.



6. Insert weight tray


7. ISTRUCTION FOR USE

Time setting

Long press PRINT key for 3 seconds to enter the TIME SETTING mode, beginning with the top row, with the flashing digit. Press PRINT key after successful change to move to next step.

EX: To input Dec 25, 2016:00am.

20 16	Enter year: Using TARE keys to decrease the value, press HOLD key to move to next digit.	
	Press PRINT key after successful entry to move to next stage.	
12.25	Enter date: Using TARE keys to decrease the value, press HOLD key to move to next digit. Press PRINT key after successful entry to move to next stage.	
08:00	Enter time: Using TARE keys to decrease the value, press HOLD key to move to next digit. Press PRINT key after successful setting.	



Display Format YYYY→MM.DD→TT:SS



ATTENTION!

For proper charging of the battery pack, connect the instrument to a power outlet for at least 8 hours

WEIGHING

Switch on the scale using key. The diagnostic scale self-check is performed and the software version is displayed.

The "0,00 kg" weight displays on the screen, now the scale is ready for weighing.

Note: If display is not at 0.0kg while switching power on. Kindly press key to zero the scale which can be used at any time to zero the scale.

HOLD FUNCTION

BABY02 is provided with the integrated hold function (determination of average value). It enables baby to be weighed accurately even they keep moving.

Once the HOLD key is pressed, the weight reading will remain on the display after the item has been removed from the scale so the reading can easily be read.

Switch 'ON' the scale using . The diagnostic self-checks is performed. The scale is ready for weighing when the "0.0 kg" displays on the screen.

Place the baby in the center of the tray gently.

Press . When the triangle is flashing on the display, the scale records the fluctuating weight values and then calculates average weight and displays the result on the screen.

Remove the baby from the scale tray. The baby's weight reading will remain on the display.

Press again to return the scale to the normal weighing mode.

HOLD function can be activated before or after putting the weight on the baby tray. But in case of weighing unstable, it is recommended to press before put baby on the baby tray.

TARA FUNCTION

Tare allows the user to zero the instrument to cancel the weight of a container/ clothes from the reading of the instrument, thus giving the true weight of the product/ person being tested.

Turn on the scale using kev.



displays on the screen and scale beeps twice, Place the object (clothes) need to be tare on the tread platform.

Press after the weight stabilizes and stable sign appears on the display.

Display will return to 000

Place the item/ baby need to be weighed (without removing the tare object). The weighing result is the Net weight of the item / baby.

To delete the saved tare value, remove the tare object from the tray and press



Using memory function (Milk difference)

This function allows automatically calculates the amount of Milk consumed

After weighting, the top row display will show the weight, for example: 3.000kg.	nold net weight
Press key to store the first weight, the second display will show "m1".	SX Loon on Model
And then press again to return to normal mode.	difference 0
When the second weight be stored, (for example: 3.180kg), press key	3.180 , s
the display will show the difference between the first and second weight.	difference
Press key to go back to normal mode, and clear the "m1" data.	<u> </u>

Note:

- This function is working only with continues on power. To cancel this function press



- By AC adapter: the scale still ON always.
- By battery: after not using of 3 minutes, the scale will OFF automatically, please keep touch by hand the
- If the difference value is over then 1000g, the display will show

8. PRINTER FUNCTIONS

Weight can be printed for records using RS232 interface cable.

After weighing simply press PRINT key to print out the results.

The format presented below is the standard format of results print-out and cannot be changed.

More informations: info@wunder.it or service@wunder.it

GROSS WEIGHT 10.00kg

TARE WEIGHT 1.00kg

NET WEIGHT 9.00kg

Set parameters of the scale interface on the connected device. It is not possible to change the scale parameters.


Baud rate: 9600 bps

Parity check: None Data length: 8 bits Stop bit: 1 bit

Handshake: RTS/CTS

Data code: ASCII





PIN	SIGNAL
2	TX
3	RX
5	GND

Serial RS232

Connector RS232

Connecting with PC

Start Hyper Terminal

Start Hyper Terminal program from clicking

Start Menu \rightarrow Programs \rightarrow Accessories \rightarrow Communication \rightarrow Hyper Terminal.

New Connection Description

Give new connection a name then click OK.

Select Your COM Port

Click Connect to select your COM port. Usually there's only one option for select. Then click OK.

Port Settings

Click Bits per second to set up rate at 9600, Data bits at 8, Parity at None, Stop bits at 1 and Flow control at Hardware. Then click OK to complete your setting.

Output Data


When press the PRINT key to output data from scale to PC or an Optional Printer.

Instruction for charging and connecting

If _____ prompt displays on the LCD, please charge the scale with exclusive adaptor.

Locate adaptor plug on the rear side of indicator.

The battery should be recharged at least every 3 months regardless if it is used or not. After a long period in storage, e.g. over 3 months, the battery should run a full cycle (charge/discharge) to allow it to restore to full capacity.



9. ERROR MESSAGE

ERROR MESSAGE	REASON	ACTION
La	Low Battery: This warning shows that the voltage of battery is too low to use	Please replace a new battery or plug the AC adaptor for operation.
Err	Overload: The total load exceeds the maximum capacity of scale	Please reduce the loading and try again.
Err.H	Counting error (too high): Indicates that the signal from the load cell/s is too height	This error is normally caused by a serious fault on the scales such as a faulty load cell or wiring. Please contact the local service representative.
Err.L	Counting error (too low): Indicates that the signal from the load cell/s is too low	This error is normally caused by a serious fault on the scales such as a faulty load cell or wiring. Please contact the local service representative.
00000	Zero count over calibration zero range +10% while power on	Please re-calibrate the instrument.
00000	Zero count under calibration zero range –10% while power on	Please re-calibrate the instrument.
Err.P	EEPROM Error: Indicates that there is a fault with the scales software	This error is normally caused by a serious fault on the scales such as a faulty load cell or wiring. Please contact the local service representative.

10. MAINTENANCE AND ASSISTANCE

For better and longer duration of the product it should receive thorough general cleaning periodically.

The instrument must be cleaned with a soft cloth moistened with water or neutral detergent, without using solvents or abrasive substances.

If the instrument remains idle for a long period, remove the batteries from the terminal.

During shipping, make sure not to subject the instrument to blows or excessive mechanical stress.

In case of repairs or assistance, contact your dealer or an authorised service centre contacting service2@wunder.it or sales@wunder.it.

11. CONFORMITY

WUNDER MODEL BABY02 ELECTRONIC SCALE SERIAL N°.....

We hereby certify that this instrument has been inspected and has successfully passed the functional test. It complies with the following standards and directives:

EN 45501 / EN60601-1-2 / EN60601-1

93/42/EEC - Medical Devices Directive

12. IDENTIFICATION LABELS



13. SCRAPPING AND WASTE DISPOSAL

If set aside for a long period, protect those parts which could be damaged due to dust build-up Scrapping

When you decide to no longer use this item, we recommend making it unusable. We also recommend making those parts which could be sources of danger harmless



EU 2002/96/EC

This product complies with the EU Directive 2002/96/EC. The symbol of the crossed-out waste bin on the appliance indicates that the product, needing to be treated separately from household waste, at the end of

its useful life must be completed in a separate collection facility for electric and electronic appliances or returned to the dealer upon purchase of a new equivalent appliance. The user is responsible for bringing the appliance to an appropriate collection structure at the end of its life. Appropriate separate collection and sending the appliance for recycling, treatment and environmentally compatible waste disposal contributes to avoid possible negative effects on the environment and health and favours the recycling of the materials the product is made of.

For more detailed information regarding available collection systems, contact your local waste disposal service or the shop where the product was purchased.

As consumers, you are obliged by law to return used or dead batteries. You may deposit old batteries at public collection spots in your town or else with any battery dealer who has placed specific collectors for this purpose. Even when scrapping electric and electronic appliances, they must be removed and deposited in specific collectors.

NOTE: The following symbols indicate the presence of harmful substances

Cadmium Mercury



ATTENTION!

Do not throw electric parts and used batteries away with household waste. Dispose of the batteries by means of your closest collection centres.

14. WARRANTY

This certificate must be kept until the warranty has expired.

It must be presented together with the invoice, tax receipt or delivery note providing the name of the dealer and date of purchase whenever a technical intervention is required. Otherwise the user will lose any warranty rights.

The warranty takes effect from the date of purchase and is valid during the entire period foreseen by the current catalogue/pricelist.

By warranty we mean the replacement or repair free of charge of parts making up the appliance which, at the discretion of the manufacturer, are deemed faulty from the origin; Wunder therefore has the faculty of repairing or replacing the item.

The warranty does not cover:

- shipping faults, damage caused by falls, carelessness or tampering
- damage caused by incapability of using the appliance and of its improper use
- damage caused by an insufficient or inadequate electrical system or alterations resulting from environmental, climatic or other types of conditions
- damage due to incorrect installation of the appliance and repairs carried out by unauthorised personnel
- Interventions at home for convenience controls or presumed defects

- Routine maintenance and that which can be considered normal wear from use
- consumables such as: power supplies, batteries, keyboards, plates, wheels, heads, rolls, load cells faulty due to blows or overloads

Service can also be refused when the appliance has been changed or transformed in any way.

In case of interventions at one's home, the customer must pay the fixed fee; if however the appliance is repaired at an authorised Wunder Service Centre, expenses and relative travel risk are borne by the user.

Wunder will not be held liable for damage of any nature caused directly or indirectly to persons, animals or objects resulting from failure to comply with all the instructions indicated in this manual or anyway resulting from improper use.

The Court of Bergamo has jurisdiction in case of any dispute

15. INSTALLATION BABY HEIGHT ROD



For the installation of digital height rod proceed as follow:

Fig.A	
Fig.B	
Fig.C	

16. INSTALL SECURITY BRACKETS

