



SC 240 MA

PORTABLE BODY COMPOSITION ANALYSER

The Worlds First Lightweight, High Capacity Medically Approved Body Composition Monitor

- Total Body Composition readings including: Weight, Body Mass Index, Fat %, Body Water %
- Additional body composition features also available when connected via software
 - Fat Mass, Fat Free Mass, Body Water Mass, Muscle Mass, Bone Mass, Visceral Fat Level, Basal Metabolic Rate, Metabolic Age
- USB Interface Connection
- Max weight capacity 200kg x 0.1kg
- Highly portable monitor, weighing only 4.7kg
- MDD Approved NAWI Class III

