Electronic wheel chair scales M500020-03









Low level platform scale with one ramp for convenient access and a foldable handrail. Two castors and one handle permit a quick transport. Powered by mains and rechargeable batteries. Allowing the use of the scale at any place.

Technical specification:

Capacity: 300kg total capacity patient

Graduation: 100 g

Dimensions: 600 x 850 x 960mm

Weight; 34kg

Power supply: By mains and rechargeable batteries

Functions: Automatic zero setting, tare, pre-tare, low battery indication,

auto witch off, hold, BMI, RS232-Interface