


NEW



Mobile wheelchair platform scale with handles and rollers

Features

- Especially suitable for weighing patients in wheelchairs because of the low-profile platform which can be approached from either side
- 1 The large steel weighing plate also makes it ideal for weighing obese patients
- Also ideal for rapid weighing of, for example, laundry carts, container trolleys, roller containers, etc. to 600 kg
- Lightweight, compact construction for mobile use
- 2 The scale can easily be transported using the two rollers and two handles and does not require much storage space
- PRE-TARE function for deduction of a known weight, for example of a wheelchair. Memory for ten PRE-TARE values
- 3 Level indicator to level the balance precisely
- Wall mount for display device, standard
- Universal external mains adapter included with delivery
- Protective working cover included with delivery

Technical data

- Large backlit LCD display, digit height 25 mm
- Dimensions weighing surface WxD 800x800 mm
- Dimensions of display device WxDxH 235x114x51 mm
- Cable length of display device approx. 5 m
- Overall dimensions WxDxH 1100x1000x80 mm
- Optional battery operation, 4x1.5 V AA not included in scope of delivery, operating time up to 60 h
- Net weight approx. 40 kg

Accessories

- Protective working cover over the display device, scope of delivery: 5 items, KERN EOB-A04BS05
- Stand to elevate display device, height of stand approx. 1000 mm, KERN EOB-A02B

STANDARD

UNIT

IP 67

BATT

MULTI

2 DAYS

1

Model	Weighing capacity	Readability
	[Max] kg	[d] kg
KERN		
MWC 600K-1	600	0,1