

DC 360 S

Dual Frequency Body Composition Analyser with Integrated Printer



NEW

Key Features

- ◊ Full body composition analysis provided in 15 seconds using clinically accurate Tanita Dual Frequency BIA Technology
- ◊ Integrated printer provides instant read out of results with topline analysis.
- ◊ Results automatically stored on the SD Card, sent to a PC or printed.
- ◊ GMon Pro Software compatible, allowing trend analysis, health risk assessments and full data management.
- ◊ Low profile platform for additional stability
- ◊ Max weight capacity 270kg with 100g accuracy

Measurements

- ◊ Body fat %
 - ◊ Fat mass kg,
 - ◊ Fat free mass kg
 - ◊ Muscle mass kg
 - ◊ Total Body Water %
 - ◊ Body mass index
 - ◊ Bone mass kg
 - ◊ Physique rating
 - ◊ Visceral fat rating
 - ◊ Basal Metabolic Rate kcal
 - ◊ Basal Metabolic Rate indicator
 - ◊ Metabolic Age
- ◊ **Print Out Analysis:**
 - Body Fat Analysis
 - Muscle Mass Indicator
 - BMR Indicator
 - Physique Rating
 - Target: BF and Weight



GMON
SOFTWARE



**WORLD No.1
BIA BRAND**

**5 YEAR
GUARANTEE**