

Model	Code	Capacity	Division
VRB-250 and VRB-300	00071I	200/300kg	100g
WHM210D	00051D	120~210cm	1mm
WH210M	00069AA	120~200cm	1mm



RB: NEW INDICATOR WEIGHT, HEIGHT, BMI

RB: MECHANICAL OR DIGITAL HEIGHT ROD



COMPACT SCALE
275x295mm

FUNCTIONS

Weight, Tare, Body Mass Index, Hold, Weight mother-child, Autooff programmable.

PRODUCT DETAILS

Certifications: ECC directive NAWI approved medical class III. ECC directive MDD Risk class Im, with measurement function.

Structure: Alluminium and metal tubular painted, ABS display, platform anti-slip.

Dimension mm: Total (w)275x(L)295x(h)1000, Platform (w)275x(L)295x(h)55, Pole (h)Ø 60x920, Display Ø(W)215x(L)180.

Display: 1° Display weight 5 digit LCD 20mm; 2° multiline display Height/BMI 3 digit LCD 12mm.

Power: External adapter and rechargeable batteries (40 hours life).

Weight: Approx 10kg

Accessories included: Adapter, Level Bubble, Adjustable feet, RS232.

Accessories not included: Electronic Height rod 100~210cm code 00051D (automatic connection height to the display); Mechanical height rod 60~200cm code 00069AA, Thermal printer WS code 01131T

Packaging: 1 pcs dimension (w)400x(L)1100x(h)210mm, Weight approx 10kg.

Warranty: 2 years.