



## BC 545N

**THE ONLY SEGMENTAL BODY COMPOSITION MONITOR IN THE WORLD WITH 10 MEASUREMENTS FOR A FAST AND DETAILED ASSESSMENT**

### Benefits

- 10 body composition measurements provided in 15 seconds
  - A unique body assessment providing essential information on the progress of the effectiveness of your training
- Muscle and Fat % measurements given for your arms, legs and trunk area displayed in easy-to-read graphics
  - Fine tune your training program and set targets for increasing muscle and decreasing fat in your arms, legs and trunk area

### Key features (compared to BC545)

- High resolution colour display including new improved easy-to-read analytical graphics
- Twice as fast measurement time - just 15 seconds
- Modern design with convenient retractable hand electrodes
- Auto recognition feature allows user to just step on for a reading
  - no buttons to press!
- Includes body mass index and visceral fat in 0.5 steps

### Measurements

- Weight
- Body Fat % (5 - 99 years)
  - Body Fat Healthy Range
  - Indicator (5 - 99 years)
- Total Body Water %
- Muscle Mass
- Physique Rating
- Bone Mineral Mass
- Basal Metabolic Rate
- Metabolic Age
- Body Mass Index
- Visceral Fat
  - Visceral Fat Healthy
  - Range Indicator
- 5 Segmental Fat Readings
- 5 Segmental Muscle Readings

### Specifications

- Auto Recognition
- Recall Function (weekly and fortnightly up to 12 months)
- Clock, Calendar
- Guest Mode
- Athlete Mode
- 5 Person Memory
- 150kg

