

# RD-545

## Connected Segmental Body Composition Monitor



Introducing the ultimate segmental body composition monitor. The RD-545 delivers the highest amount of accurate readings with dual frequency - which can only be found in professional medical devices.

### Benefits

- Highly accurate through dual frequency segmental reactance/resistance technology
- More than 20 body composition measurements
- 50 gram graduation
- Integrated Bluetooth technology
- Connected through “My Tanita App” for Android and I-phone
- Recognition mode, just step on
- 4 Memories and color coded indicators
- 5 years warranty



## Measurements

Your whole body will be analysed in detail using dual frequency resistance/reactance technology on segmental level. Arms, legs and Trunk are measured separately.

Changes and individual differences for all readings can now be analysed with medical accuracy.

- ◊ Weight in kg
- ◊ Body fat %
- ◊ Muscle mass in kg
- ◊ Body mass index (BMI)
- ◊ **NEW!** Muscle quality
- ◊ **NEW!** Muscle score
- ◊ Physique rating
- ◊ Bone mass in kg
- ◊ Total body water %
- ◊ Visceral fat
- ◊ Basal metabolic rate
- ◊ Metabolic age
- ◊ Segmental measurement of left and right arm, fat and muscle mass
- ◊ Segmental measurement of left and right leg, fat and muscle mass



World's first dual frequency segmental body composition monitor with 8 electrodes measurement technology connected with My Tanita app.

## Technology

- ◊ BIA, dual frequency 50kHz, 6.25 kHz, 100uA

## Features

- ◊ Integrated Bluetooth technology
- ◊ 5 year warranty
- ◊ High accuracy
- ◊ 4 x A4 batteries included
- ◊ 200 kg

## EAN code

- ◊ 4904785817029

