MB120 with SmartGuide™



SmartGuide™ get started quickly


This revolutionary feature automatically analyzes a sample and creates a method,

thereby removing the guesswork from the initial setup

- With or without a reference moisture value (from the oven method or other reliable source)
- Possibility to save the created method directly to the MB's memory

User Management

 MB120 provides operational control and aids compliance by restricting access so that unauthorized users may use the device but not change settings



Halogen heating

 MB is capable of producing fast and accurate results due to features such as a halogen heating system, customizable shut-off criteria and four built-in drying profiles

Method and results memory

 MB120 allows for storing up to 100 methods and 1000 measurement results. Stored results can be exported and methods can be shared across multiple devices

