



PORTABLE PERSONAL SCALE VB3-200

The Scale for health authorities, athletic games, wrestling, boxing, etc.

FEATURES

- Capacity 200 kg (440 lb). Graduation 100 g (0,2 lb)
- High quality and excellent accuracy
- Robust performance in aluminium
- Alkaline batteries and adapter included
- Kg/lb switch
- Auto power off
- Swedish design
- Nylon case included

OPTIONS

- RS232
- Blue tooth
- Wlan
- USB
- Printer MP-20

SPECIFIKATIONER

Supply:	6 x AA alkaline batteries or included external 9VDC adaptor
Dimension, weight:	260 x 345 x 55 mm, 4,9 kg
Functions:	On, Net/Gross, Off/Zero, Tare, Auto-power off
Data interface:	RS232
Temperature Range	-10 to +40°C (Approved models 0 to +30°C)

MODEL	CAPACITY/GRADUATION	EC-TYPE APPROVAL	ANRN
VB3-200-EC	Max 200 kg/ 100g	Yes	FI01.1.03. Revision 1 according to 90/384/EEC. VB2-200-EC Medical Directive 93/42/EEC CE Accuracy class: III
VB3-200-50	Max 200 kg/ 50 g	No	"Sport Scale" for all applications without trade and medical use.