


DATASHEET WBA300



New indicator Weight, BMI, Fat Mass



BODY ANALYSIS SCALE



Model	Code	Capacity	Division [e]	% Fat Mass
WBAcolonna	00052D	300kg	100g	0,5
Stampante WS	01131T	—	—	
Stat.HM210D	00051C	120~210cm	1mm	
Stat.Mecc.	00069AA	60~200cm	1mm	

Ecc directive NAWI Approved medical class III, for medical use

Structure Aluminium and metal tubular painted, ABS display, 4 stainless electrode, 2 wheels (easy to carry)

Functions Weight, Tare, BMI, Fat mass %-kg, Free Mass kg; Total Body Water kg; Basl Metabolic

Display ABS indicator with DOUBLE LCD triple reading 1° Weight 5 digit LCD 20mm only Weight; 2° multiline display Height/BMI and Fat/Age 3 digit LCD 12mm

Power External adapter rechargeable batteries (50 hours life)

Size mm	Total (W)350x(L)600x(H)920	Platform (W)350x(L)450x(H)85	Display Ø(W) 215x(L)180	Pole (H)920
---------	-------------------------------	---------------------------------	----------------------------	----------------

Weight Approx 10kg

Packaging 1 pcs dimensions (w)400x(L)1100x(h)210mm , Weight approx 22kg

Accessories included RS232, Adapter, 2 Wheels to move, Level Bouble, Adjustable Feet

Accessories not included Electronic Height rod 120~210cm code 00051C (automatic connection height to the display)
Mechanical height rod 60~200cm code 00069AA,
Thermal printer WS code 01131T

Warranty 2 years