

PRODUCT INFORMATION



Class III, approved Scales
Electronic Wheelchair and Multi-purpose Scales



Wheelchair Scale

- Large LCD, 25 mm
- Automatic zero setting and function test
- Tare and hold function
- Pre tare setting function, incl. 10 memories
- BMI (Body Mass Index) function
- Power supply: mains and rechargeable batteries
- Print function and PC connection via serial interface RS 232
- Can be stored vertically


Options:



Easy manoeuvrable on two castors

Options:

2nd Ramp, side boundaries & foldable column



Technical Specifications

Capacity	300 kg
Graduation	100 g
Dimensions	915 x 44 x 800 mm
Weight	approx. 22 kg
Platform size	800 x 800 mm



M500020