




More information on the website
radwag.com/en/info,w1,Q01

Fingerprint Reader



Overview

An accessory for the ELLIPSIS laboratory balances, designed to read fingerprints in order to log into the balance. Simply touch the reader to log in. The fingerprint reader is connected to the ELLIPSIS balance via the USB port.



Compatible with

XA 5.M.A.P Microbalance
 Microbalance for Pipette Calibration MYA 5Y.P
 PM 5Y Precision Balance
 Microbalance for Filters MYA 5Y.F
 XA 5.M Microbalance

XA 5Y Analytical Balance
 UYA 5Y Ultra-Microbalance
 MYA 5Y Microbalance
 Ultra-Microbalance for Filters UYA 5Y.F