

GMON Professional

TANITA'S GMON PROFESSIONAL SOFTWARE IS DESIGNED TO LINK EASILY VIA RS-232C OR USB TO ANY SCALE IN THE TANITA PROFESSIONAL RANGE.

The health monitoring software can be used to chart progress and trends in composition measurements over time and includes the facility to produce professional reports and printouts. Client information can be stored with analysis data to give personal trainers, healthcare professionals and serious athletes accurate, in-depth information to help them devise tailored training, fitness and weight-loss programmes.

GMON Professional is ideal for:

- Health screening
- Weight management
- Fitness assesment

System Requirements

- Windows 2000 / XP / Vista / Windows 7 / Windows 8
- Monitor minimum size: 1024 x 768 x 24 Bit

Product Details

- Direct data transfer to and from all professional Tanita scales
- Direct data transfer to and from the Activity Monitor AM-180
- Manual input for all scales
- Manual input and analysis of other body and vital values
- Calculation of personal health risk
- Target Setting
- Professional reports and printouts
- Database query, extract whatever data you wish and save it in XLS, CSV, TXT, ... format
- MDD Class I Certified

