

## TP2100 Thermal Printer

Connect thermal printer to compatible Charder scales and print out measurement data from display panels including weight, height, BMI, tare and prietare weight with one button. Printer is connected to device via USB cable or Wireless transmission (TP2110 only).

## Specifications

Overall Dimensions	I20(₩) x 80(D) x 50(H) mm
	Number of dots: 384 dots/line Resolution: 8/dots/mm
Printing	Paper: 57.5+0.5mm(W)x@30mm
	max
	Printing speed: 75 mm/s
Device Weght	0.2 kg
Data Transfer	RS232 to USB



Power Supply I 2V Power Adapter

Compatible devices DP2701, DP3400, DP3710, DP4400, DP4600, DP4600, DP4800D