

## WB 380

### DIGITAL SCALE WITH BMI FUNCTION

The WB380 is a high capacity precision digital scale featuring a BMI function. The scale is available in three formats: portable, with a column mounted display or with an integrated height rod.

#### At a glance

- High weight capacity 300kg
- Step on activation
- Low level, stable platform for elderly or obese patients
- Swivel easy-to-read LCD display allowing confidential readings

#### Measurements

- Height
- Weight
- BMI



Model		WB-380S	WB-380P	WB-380H
		Remote Display	Column Mounted	With Height Rod
Weight Measurement	Maximum Capacity / minimum Graduation	300kg / 0.1kg		
Height Measurement	Measurement System			Mechanical Height Rod
	Range of Height			64cm – 214cm
Input Items	Height	61cm- 250cm 1cm increments		
Output Items	Display	Weight	300kg / 0.1kg	
		Height	61cm- 250cm 1cm increments	
		BMI	0.1 increments	
Product Size	Overall	240 x 139 x 123 mm (Indicator)	395 x 595 x 1196 mm	395 x 552 x 1425 mm
	Platform	395 x 390 x 59 mm		
Product Weight		Total 7.1kg	Total 10.0kg	Total 11.2kg
Output Data Interface		RS-232C (D sub 9-pin Female connector) USB (B-type)		
Power Source		AC adapter (included) Centre Minus DC 6V 200mA (LR6 - AA Alkaline Battery x 4)		
Battery Life		Approximately 100hours of continuous use		