





### The triple beam that truly has it all!

Accurate, dependable and stylish; the OHAUS Triple Pro<sup>®</sup> mechanical balance offers it all: The superior quality and durability you've come to expect from OHAUS, plus all the extras to satisfy every classroom activity requiring a balance. Just like the original Triple Beam balance from OHAUS, the Triple Pro balance gives you the confidence to know you have the accuracy and convenience of a toploading balance while allowing your students to gain hands-on experience through manipulation of the poises, attachment weights and the rod & clamp assembly (all included). The OHAUS Triple Pro balance is the perfect learning tool to help students visualize mass measurement principles and to emphasize a hands-on approach to learning.

### **Standard Features Include:**

- Unique, easy grip carrying handle is ergonomically designed and perfectly balanced to allow for effortless movement around the classroom
- Two security features are included to prevent theft; a security loop and a Kensington-type lock (cables sold separately)
- Includes a convenient on-board storage compartment for the included attachment weight set
- A built-in transportation/storage lock helps to protect the longevity by minimizing wear on the knives and bearings, and reducing the possibility of bending the damper vane
- Includes rod & clamp assembly for below-balance weighing (specific gravity determination)
- Backed by an incredible 10 year limited warranty that only OHAUS can provide

# **Triple Pro<sup>®</sup>** Mechanical Balance

## **Standard Features:**

Readability

Calibrations:


Front beam

Rear beam

Accessories

Middle beam

Stainless steel plate platform



0.1g

10 g x 0.1 g

500 g x 100 g

100 g x 10 g 5.8 x 5.8 inches / 14.73 x 14.73 cm

**Ohaus Number** 

www.ohaus.com

With offices throughout Europe, Asia, and Latin America

#### \* ISO 9001:2000

© Copyright Ohaus Corporation

80774440



