



BC 601

SEGMENTAL BODY COMPOSITION MONITOR WITH SD-CARD

Exclusive Features

- · Provides segmental body composition data for each arm, leg and trunk area
- Includes an SD card to transfer all results to PC
- Free software enables user to view and manage data
- · Hand held, pull-up console features clear LCD display, and set-up buttons

Specifications

Weight Capacity: 150kg

· Graduation: 100g

Display:45 x 125mm

· Memories: 4

· Batteries: Included

Guarantee: 3 Years

 Modern design with convenient retractable hand electrodes

Measurements

- Weight
- Body Fat % (5 99 years)
 - Body Fat Healthy Range Indicator (5 - 99 years)
- Total Body Water %
- Muscle Mass
- Physique Rating
- Bone Mineral Mass
- · Basal Metabolic Rate
- Metabolic Age
- Body Mass Index
- Visceral Fat
 - Visceral Fat Healthy Range Indicator
- 5 Segmental Fat Readings
- 5 Segmental Muscle Readings

