

# VRH-300 HANDRAIL OBESITY SCALE



<i>Ecc directive NAWI</i>	<i>Approved medical class III, for medical use</i>			
<i>Ecc directive MDD</i>	<i>Risk class I, with measurement function</i>			
<i>Structure</i>	<i>Alluminium and metal tubular painted, ABS display, platform anti-slip, 2 wheels (easy to carry)</i>			
<i>Functions</i>	<i>Weight, Tare, Body Mass Index, Hold, Weight mother/child, Autooff programmable, PRINT: WEIGHT, TARE, HEIGHT, BMI (BODY MASS INDEX)</i>			
<i>Display</i>	<i>1° Weight 5 digit LCD 20mm; 2° multiline display Height/BMI 3 digit LCD 12mm</i>			
<i>Power</i>	<i>Adapter rechargeable batteries (50 hours life)</i>			
<b>Dimensioni mm</b>	<i>Total</i>	<i>Platform</i>	<i>Display</i>	<i>Colonna</i>
<i>Size mm</i>	<i>(W)550x(L)550x(H)1160</i>	<i>(W)550x(L)550x(H)85</i>	<i>Ø(W)215x(L)180</i>	<i>(H)920</i>
<i>Weight</i>	<i>Model RH approx 20kg; Model RO approx 14kg</i>			
<i>Packaging</i>	<i>1 pcs dimensions (w)400x(L)1100x(h)210mm , Weight approx 22kg</i>			
<i>Accessories included</i>	<i>RS232, Adapter, 2 Wheels to move, Level Boule, Adjustable Feet</i>			
<i>Accessories not included</i>	<i>Electronic Height rod 120~210cm code 00051C (automatic connection height to the display ) Mechanical height rod 60~200cm code 00069AA, Thermal printer WS code 01131T</i>			
<i>Warranty</i>	<i>2 years</i>			

