

# PRODUCT INFORMATION

# ADE

Class III, approved Scales  
Electronic Bed Scales



**M601020**

## Bed Scale

- Large LCD display, 25 mm
- Automatic zero setting and function test
- Tare and hold function
- Pre tare setting function incl. 10 memories
- BMI (Body Mass Index) function
- Power supply: mains and rechargeable batteries
- Plug connection to load cell extension lead
- Print function and PC connection via serial interface RS 232

Two independent weighbridges with  
built-in-ramps for easiest positioning of bed



### Technical Specifications

Capacity	500 kg
Max weight of patient	300 kg
Graduation	100 g
Dimensions	1,000 x 90 x 450 mm
Weight	approx. 29 kg