


PRODUCT INFORMATION

ADE

Class III, approved Scales
Electronic Bed Scales



M600020

Bed Scale

- Large LCD, 25 mm
- Automatic zero setting and function test
- Tare and hold function
- Pre tare setting function, incl. 10 memor
- BMI (Body Mass Index) function
- Power supply: mains and rechargeable batteries
- Print function and PC connection via serial interface RS 232
- Incl. mobile stand
- Safe adapter holder

Four weighing shoes with
integrated ramps

Technical Specifications

Capacity	500 kg
Max weight of patient	300 kg
Graduation	100 g
Dimensions	600 x 920 x 550 mm
Weight	approx. 33 kg